

supplier diversity portfolio

A collection of quality products sourced from diverse suppliers and backed by the reliable support and service you expect from Owens & Minor

2010



MEDICHOICE[®]
from Owens & Minor

Your choice for quality, value and unmatched service



MIX
Paper from
responsible sources
FSC® C051259

The manufacturers of the products contained in this catalog are the sources of all information herein relating to the products, including, but not limited to, product content and warnings. Owens & Minor makes no representations or warranties expressed or implied regarding the product information contained herein. If further information about any product is desired, please contact an Owens & Minor representative, who will assist in obtaining the information from the manufacturer.

2010 supplier diversity portfolio

Advancing diverse sourcing through partnerships

Owens & Minor has had a long-standing commitment to advancing supplier diversity in healthcare. We have a recognized track record of developing diverse relationships and of championing processes for collaboration that build infrastructure and ethical business practices within a diverse group of woman-, minority-, and veteran-owned manufacturing and service organizations.

We believe our investment in supplier diversity and supplier development promotes a superior supply chain for our provider customers. We also know that a thriving community of diverse suppliers generates healthy competition and innovative delivery methods that can help drive cost out of the supply chain.

Supplier diversity and MediChoice®

As Owens & Minor's MediChoice® line has grown, our commitment to supplier diversity in our sourcing relationships has grown with it. In fact, today over 30% of MediChoice® product sales are sourced from woman-, minority- or veteran-owned businesses. Diversity has become a strategic advantage for MediChoice®. By choosing MediChoice, hospitals not only achieve significant savings on their medical-surgical supply purchases, but also they can advance their own supplier diversity goals.



“Commitment to supplier diversity is an integral part of Owens & Minor’s culture. We recognize that fostering relationships with minority-, woman-, and veteran-owned businesses drives innovation, deepens our understanding of our customers, fuels mutual growth, and strengthens the communities we serve.”

—Gavin Jeffs, Operating Vice President, Supplier Relations

admissions

PATIENT SUPPLIES.....10
 PERSONAL CARE10
 PILLOWS12

anesthesia

Berman Airways 14

apparel—patient

Mesh Briefs15
 SLIPPERS15
 TELEMETRY POUCHES17

apparel—staff

FACE MASKS18
 ISOLATION GOWNS18
 DISPOSABLE COVERALLS 21
 DISPOSABLE LAB COATS 22

central supply

Enzyme Gel Spray Pre-Cleaner 23
 Enzyme Cleaner Pre-Soak..... 23
 Instrument Lubricants 23
 Neutral pH Instrument Cleaners 23
 INSTRUMENTS 24

diagnostics

Lubricating Jelly 25
 HAND-HELD SPHYGMOMANOMETERS .. 25
 STETHOSCOPES 27
 DIGITAL THERMOMETERS 31
 EKG Paper 32

equipment

BLANKET WARMERS 33

kits & trays

IRRIGATION TRAYS 35
 Incision and Drainage Trays 36
 Central Line Dressing Trays 36
 IV START KITS 38
 Suture Trays 42
 Skin Staple Remover Trays 43
 Suture Removal Trays 44
 Tracheostomy Trays..... 45

maternity and baby care

Fetal Monitoring Belts..... 46
 Baby Bonnets 46
 Baby Lotion..... 47
 Baby Shampoo 47
 Perineal Cold Packs 48
 PERINEAL WASHES..... 48

O.R. products

Biopsy Punches.....49
 Disposable Safety Scalpels.....49
 Esmark Bandages.....50
 SKIN MARKERS50
 Suction Tubing.....51
 Tubing Connectors.....51
 Yankauers.....52

patient care

Alcohol and Peroxide.....53
 Instant Hand Sanitizers53
 Instant Cold Packs.....54
 Instant Warm Packs54
 Reusable Hot or Cold Gel Packs55
 Adult Wipes55
 Washcloths56
 Odor Eliminators.....56
 Oral Swabs57
 Tourniquets57

respiratory

Suction Catheters58
 Suction Catheter Kits.....58
 Suction Catheter Trays59

urology

Irrigation Syringes.....60
 Pediatric Urine Collectors60
 Foley Insertion Trays61
 Urethral Catheterization Trays61

wound care

Cotton Balls.....62
 Cohesive Bandages.....62
 ELASTIC BANDAGES.....63
 Tubular Elastic Support Bandages.....64
 Tubular Net Bandages.....65
 Paste Bandages66
 Compression Bandage System.....66

product index

ALPHABETICAL INDEX 67

Look for the following symbols throughout the catalog



Indicates Latex-Free



Indicates Diverse Supplier



Indicates Ethylene Oxide Sterilized



Indicates Gamma Radiation Sterilized



Indicates Aseptic Processing Sterilized

Amsino International, Inc.

Amsino International, Inc., established in 1993 by Dr. Richard Y. Lee, is a leading developer and manufacturer of disposable medical devices. The California-based, ISO-certified company produces IV/infusion, surgical, urological, enteral feeding/irrigation and patient care products, which are distributed globally. Amsino, a certified small and minority-owned business, also provides full-service contract manufacturing and sourcing services.

Two thousand employees strong, the company has in-house engineering and R&D teams. Its state-of-the-art manufacturing equipment produces a broad range of injection molded products. In addition, the company has radio frequency sealing, ultrasonic welding, aseptic filling, EtO and steam sterilization, and form-fill-seal packaging capabilities. Amsino products are assembled in medical-class clean rooms.

Amsino's relationship with Owens & Minor began in 2004. Today, Amsino supplies respiratory, surgical, personal protective wear and urology products to the MediChoice product portfolio.

The company's mission is to provide complete customer satisfaction in all areas of operations: product design and development, engineering, technical support, manufacturing, quality control and inspection, shipping and logistics, sales and marketing, and customer service.

AVCOR Health Care Products, Inc.

AVCOR Health Care Products, Inc. is a veteran-owned business enterprise, founded in 1984 by Fred M. Trainor, an Army veteran with over 45 years of experience serving healthcare. Mr. Trainor founded the Hospitex Division of American Hospital Supply Corporation and TecnoL, Inc. which is now a division of Kimberly Clark.

AVCOR's initial product segment was Compression Dressings featuring the Honeycomb Elastic Bandage with a knitted design and Velcro® closures to replace woven

products with metal clip closures. The Honeycomb product soon became the premium elastic bandage market leader and it remains a leader today.

Over the years, additional innovative types of compression bandages have been added to AVCOR's line, including cohesive and economy knits and latex-free products. AVCOR introduced the first latex-free products in wound dressings and blood collection products with a unique lightly textured blue tourniquet in the late 1980s. AVCOR has sought to deliver product offerings that provide maximum nursing convenience, patient safety and therapeutic use.

AVCOR partners with Owens & Minor to provide MediChoice cohesive bandages and flat tourniquets.

Care Line Industries, Inc.

Care Line Industries, Inc. is a privately held veteran-owned, small business. The company has been a leader in providing high-quality patient care items in both bulk and kit form for over 40 years. Care Line has the experience and supply chain to meet the value-driven needs of its customers. The company started the commercial sale of non-sterile patient care kits in 1965 as a Promotion Service Company. In 1983, Care Line began manufacturing disposable and reusable bed pillows. Care Line now provides a complete selection of these product types through four manufacturing facilities.

Care Line knows the role that service plays in its success. Its customer service department is second to none and its goal is to work with customers to understand their needs and to help select products that will make the customer successful. Patient-comfort is the cornerstone of Care Line's decision making. One of Care Line's mottos is "small enough to care and large enough to serve."

Since 2004, Care Line has provided disposable, reusable and limited reusable pillows to Owens & Minor's MediChoice product line.

Cellucap

Cellucap Manufacturing Company is a woman-owned business enterprise that manufactures a full range of disposable protective apparel for the healthcare, food service, safety, industrial and food processing industries. The family-run business combines its company-owned, U.S. manufacturing capabilities with a strong import program to provide customers worldwide with a consistent, uninterrupted supply of high-quality products at competitive prices.

In 1998, Cellucap acquired Melco, a leading manufacturer of disposable garments. This acquisition strengthened Cellucap's product line, manufacturing capabilities, long-term supply chain relationships and distribution outlets nationwide.

The company provides products ensuring safety and also emphasizes safety in its manufacturing processes. For instance, its bouffant caps remain untouched during the entire production process. Cellucap works one-on-one with customers to select the proper fabric for safety gowns to provide the appropriate level of personal safety in various patient care environments.

Cellucap provides disposable staff apparel for Owens & Minor's MediChoice brand products.

Century Hosiery, Inc.

Century Hosiery, Inc. was founded by Kathy and Malcolm Martin in 1989 and is a certified woman-owned business. Century provides traction footwear for the healthcare and retail markets, manufacturing its products exclusively in the U.S. and applying strict quality control measures throughout the production process. The company's North Carolina headquarters houses knitting, dyeing, sewing, printing, warehousing and shipping operations all under one roof.

Originally, Century Hosiery outsourced knitting and finishing, but in 1995 it expanded its operations, doubled its staff to 70 and expanded its customer base. Today, Century has 140 employees and is fast becoming a leading U.S. manufacturer of medical footwear known for its reasonably priced, high-quality products and its customer service.

A partner of Owens & Minor since 2003, Century provides a full line-up of patient care slippers for the MediChoice product portfolio. In 2009, Century Hosiery was selected as Owens & Minor's MediChoice Supplier of the Year.

ColdStar International, Inc.

ColdStar International, Inc. is a woman-owned business enterprise that has provided a quality line of hot and cold therapy products for over 47 years. In 1963, Verlam Industries pioneered a new technology in cold therapy. In 1977, Sam Francis, CEO of ColdStar International purchased the assets of Verlam Industries and incorporated it as Jack Frost Laboratories, Inc. Francis insulated the first reusable instant cold pack to offer the medical community a solution for "Pennies per Application."

In the early 1980s, Sam began to implement improvements to the products and automate the production of hot and cold therapy products. Jack Frost became the first automated instant cold and hot pack company in the world. ColdStar International's enduring commitment to quality and innovation is its foundation.

Today, ColdStar International, Inc. produces a comprehensive line of the most innovative hot and cold therapy products. All products are American-made—manufactured in Northeast Tennessee. ColdStar supplies hot and cold therapy products to the MediChoice product portfolio.

Hydrox Laboratories

Hydrox Laboratories, a leading supplier to the healthcare and beauty supply industries, is a certified minority-owned business enterprise. The Illinois-based company's products are distributed to more than 5,000 hospitals and 3,500 retail establishments nationwide.

Founded in 1913, Hydrox is among the nation's oldest continuously running manufacturers. It is the U.S. patent originator of hydrogen peroxide solution, which it produced exclusively for 70 years. To serve a broader market, Hydrox

expanded its product line in the early 1980s. Today, its products include mouthwashes, deodorants, hydrogen peroxide solutions, isopropyl rubbing alcohol, lotions and moisturizers, shampoos, soaps and hand sanitizers.

The company's 80,000 square feet of production and warehouse space is custom-designed to meet world-class standards for flexible manufacturing and fast response times for which Hydrox has won numerous supplier awards. Its leading-edge technology helps ensure reliable product purity and safety.

Hydrox supplies hydrogen peroxide solutions, isopropyl rubbing alcohol, mouthwash, lotion, shampoos and hand sanitizers to the MediChoice product portfolio.

Innovate Medical, LLC

Innovate Medical, LLC, is a small, woman-owned business that was founded in 2001 by Susan F. Johnston, a Tennessee patent attorney with a background as a pharmaceutical and polymers chemist. Having served as legal counsel for medical device inventors and businesses for several years, she acquired a group of medical device patents and formed a small working team to focus on product innovation and developing customer solutions and market strategies.

Since 2002, Innovate has been the nation's highest volume contract manufacturer of latex and latex-free rubber Esmark pressure bandages. Related medical products include Texstrip Elastomer bandages, which consistently rank high in customer satisfaction and reliability, Inno-Vascular brand latex-free rubber tourniquets, and exercise bands. Innovate's product development efforts include specialty niche disposable patient care devices and kits.

Since 2006, it has provided reeled latex-free tourniquets to the MediChoice line.

International Medsurg Connection

Founded in 1996 by husband-wife team, Dr. Omi Bhati and Nalini, International Medsurg Connection (IMC) is certified as a small, minority-owned business. The company offers a product portfolio that aims to help healthcare providers manage economic constraints while delivering high quality patient care.

The company began as a cotton textile supplier of OR towels, lap sponges, infant receiving blankets and gauze dressings. With that focus, IMC built strong relationships with custom tray manufacturers and medical-surgical supply distributors. Since its beginnings, IMC's product portfolio has grown significantly. The company now offers items in more than 11 distinct product segments and has annual sales exceeding \$20 million.

IMC's core values center on quality and credibility in its product offerings. Supporting those values is a management team with over 50 years of collective experience in the medical field and a global presence with infrastructure in the United States and abroad.

IMC supplies patient care products to the MediChoice product portfolio.

Kerma Medical Products, Inc.

Kerma Medical Products, Inc. was founded in 1991 by wife-husband team, Joyce Kershaw-Reubel and Earl Reubel, as a minority and woman-owned business. In 2006, Joe Reubel became company president, making Kerma a minority and veteran-owned business.

Virginia-based Kerma originated as a manufacturer and supplier to the United States government. When Congress cancelled all federal direct supplier contracts several years later, Kerma refocused its business on non-government affiliated acute-care market. Today, the company has contracts with major group purchasing organizations (GPOs), it supplies products to over 60% of the nation's acute care hospitals, and it has regained some of the government-sector business it initially served.

In 1995, Kerma entered a Department of Defense Mentor-Protégé relationship with Owens & Minor. In 1999, the partnership earned a Nunn-Perry Award, granted by the Department of Defense (DOD), Office of Small Business Programs for outstanding mentor-protégé teams. By 2006, Kerma joined Owens & Minor's Top 50 Supplier Listing and earned Owens & Minor's Supplier of the Year Award.

Kerma manufactures soft goods, including burn dressing pads, fetal monitor belts, disposable washcloths and baby bonnets. Kerma also assembles and/or packages elastic bandages, telemetry pouches, disposable face shields, digital thermometers, stethoscopes and blood pressure cuffs. Many of these products are part of Owens & Minor's MediChoice private-label offering.

Kerma believes product selection, quality exceeding customer expectations, competitive pricing, timely delivery and excellent customer service are central contributors to its ongoing success.

MAC Medical Supply Company

MAC Medical Supply Company, a medical-surgical supply manufacturer and distributor with a national customer base, is a woman-owned business certified by the Women's Business Enterprise National Council.

The company started as provider of technically advanced products and services for the cardiology and electrophysiology fields and continuously expanded its product line in the medical-surgical arena with a strengthened focus on medical supplies. In 2000, entrepreneur Millie Maddocks acquired and incorporated the supplies division that today is MAC Medical Supply Company.

MAC Medical's product line includes a broad selection of medical and surgical supplies such as blood pressure devices and accessories, chemotherapy protection and safety products, defibrillator supplies, disposable electrodes, EKG cables and respiratory products.

In 2003, MAC Medical became a supplier to Owens & Minor, providing medical imaging and recording chart papers under the preferred supplier program. In 2007, MAC Medical became Owens & Minor's MediChoice private-label

supplier for medical recording chart papers and enzymatic instrument cleaners. In that same year, Owens & Minor and MAC Medical formed a Mentor-Protégé relationship under the Department of Defense Small Business Office, signifying a new stage of growth for MAC. This Mentor-Protégé relationship with Owens & Minor continues today.

In 2009, Ms. Maddocks was honored to be selected as a participant in the *Inc. Magazine*-William J. Clinton Foundation Partnership's 2009-2010 Entrepreneur Mentoring program.

MAC Medical holds GPO contracts with HealthTrust Purchasing Group (HPG), Consorta and the Department of Defense (DAPA). The company also holds FSS contracts with the Veterans Administration and is an authorized distributor for MedAssets.

Medical Resources

Medical Resources was founded in 1987 by Suzi Reichenbach and is a certified woman-owned business enterprise. Medical Resources is dedicated to providing quality medical equipment and supplies, technical expertise, superior service and cost effective pricing to healthcare professionals, ambulatory surgery centers and hospitals. Getting its start as an office furniture dealer, Medical Resources soon expanded its offering to include floor and wall coverings, and cabinetry and today its line also includes medical equipment and supplies.

Medical Resources manufactures stainless steel warming cabinets, instrument tables, scrub sinks, case carts, endoscope drying cabinets, storage cabinets, case carts and custom cabinetry. Medical Resources provides interior design services, space planning, standardization programs, seating, recliners, desks, cabinetry, tables, bariatric furniture, lab cabinetry, panel systems, filing, shelving, storage solutions, artwork and accessories. Medical Resources also provides automated, computerized remote temperature and humidity monitoring solutions designed to monitor, record, and alert pharmacy, dietary, SPD, blood bank, and laboratory departments as well as offsite facilities.

Medical Resources provides stainless steel warming cabinets for Owens & Minor's MediChoice brand products.

Medikmark, Inc.

Medikmark, Inc. was founded in 1987 and is a veteran-owned, small business. Medikmark is a manufacturer of sterile and non-sterile nursing care kits and medical contract packaging services. Medikmark is a leading supplier of standard and custom medical procedural kits and trays offering a complete line of products ranging from standard and custom dressing change, IV start and laceration kits and trays to IV sets, dialysis and staff protection kits, trach care trays, instrument trays and sterile solutions for a wide variety of healthcare uses.

Medikmark's FDA registered location in Waukegan, IL, contains a certified Class 100,000 clean room for procedural kit and tray production and a certified Class 10,000 clean room for sterile solution production. Its versatile production capabilities and flexible staffing enables it to respond to customers of all sizes—with small, medium and large production requirements.

Medikmark serves as a supplier of minor and custom procedure kits and trays for the MediChoice product line. MediChoice kits from Medikmark can be customized to a customer's exact specifications including contents, quantities, and brand name of components.

Medi-Tech International Corp.

Medi-Tech International Corporation has been producing a high-quality line of wound management products for hospitals, extended care facilities and physician offices since 1971. Its product development focuses on providing a more comfortable wound management experience for the patient, while also saving time and money for the healthcare provider. The company's customer base is international, yet it maintains its manufacturing operations in the United States.

In 2007, the New York-based company changed its diversity supplier status to Small Disadvantaged Business when it began functioning as a woman-owned company.

Medi-Tech invented the original tubular elastic dressing retainer net, Spandage™, and later followed suit with the MT Spandage™ line. In addition to these innovative dressing securement products, Medi-Tech also offers skin therapy washes and lotions, wound care/burn dressings, mesh briefs and a patient safety product.

Medi-Tech has been an Owens & Minor partner since 2001, and manufactures MediChoice tubular net bandages and MediChoice disposable mesh briefs.

OP-marks, Inc.

OP-marks, Inc., is a woman-owned and operated business specializing in surgical site marking products, specifically designed to help surgeons and institutions guard against wrong-site surgeries. OP-marks products are used throughout the United States and Canada.

Since its founding in 2004, OP-marks has been an innovator in the field of surgical site marking. OP-marks was the first to introduce the mini skin marker—now the industry standard site marking tool. It is the only company that supplies these markers in proof-of-use packaging to guard against cross-contamination of resistant pathogens. OP-marks was also the first company to introduce temporary skin laminates as site marking devices. These also have become the industry standards for site marking Ophthalmology and ENT patients.

In collaboration with the Georgia Institute of Technology and Research, OP-marks developed the first truly "prep-proof" ink for site marking patients. This revolutionary hybrid ink ensures that the surgical site remains clearly marked after any form of surgical skin prep.

Trinity Sterile, Inc.

Trinity Sterile Inc. acquired the healthcare kit-packing capabilities of Trinity Laboratories Inc. in 2004 to form its current company specializing in both sterile and non-sterile kits and related medical supplies. The Maryland-based company is a certified minority-owned business enterprise, with 70 percent of its employees having a minority status.

Since its acquisition, Trinity's sales have nearly quadrupled, and employment has doubled. Combined, company management has more than 100 years experience in production and distribution.

Trinity Sterile's product line includes sterile convenience kits and accessories, ranging from catheter insertion kits and tracheostomy care kits to a wide range of IV start kits. Trinity also carries non-sterile items such as tourniquets, ear cures and floor-grade instruments.

Trinity sterilizes all products in-house, to ensure both strict adherence to its quality standards throughout the manufacturing process and quick response times. Within its 45,000 square-foot facility is a Class 10,000 clean room where employees assemble all products.

Trinity Sterile supplies standard and custom kits and trays to the MediChoice product line.



PATIENT SUPPLIES
Denture Cups

- Teal bowl with opaque white lid
- Can write patient information on lid
- Convenient 8-oz. size
- 25 per pack; 10 packs per case
- Durable polypropylene construction



Item No.	Description	Quantity
33660DC5000	Denture Cup, Opaque, Teal, 8 oz. w/ Lid	1 / EA, 250 / CA



PERSONAL CARE
Hand & Body Lotion

- Available in cocoa butter fragrance and fragrance free
- Contains aloe vera and Vitamin E to promote skin health
- Quick-absorbing, non-greasy formula
- Emollients and moisturizers help replenish natural skin oils
- Convenient flip-top cap for easy dispensing
- Extra mild, pH balanced
- Made in USA



Item No.	Description	Quantity
3306PC4020	Hand and Body Lotion, 2 oz., Scented	1 / EA, 144 / CA
3306PC102FF	Hand and Body Lotion, 2 oz., Fragrance Free	1 / EA, 144 / CA
3306PC4035	Hand and Body Lotion, 4 oz., Scented	1 / EA, 96 / CA
3306PC104FF	Hand and Body Lotion, 4 oz., Fragrance Free	1 / EA, 96 / CA



PERSONAL CARE
No-Rinse Shampoo & Body Bath

- Gentle shampoo cleans without soap and water
- Can be used on hair or skin
- Convenient flip-top cap for easy dispensing
- Gentle cleansing, less irritating than soap
- pH balanced
- Made in USA



Item No.	Description	Quantity
3306PC9020	Shampoo, No-Rinse, 2 oz. w/ Flip-top cap	1 / EA, 144 / CA
3306PC9080	Shampoo, No-Rinse, 8 oz. w/ Flip-top cap	1 / EA, 48 / CA



PERSONAL CARE
Roll-On Deodorant

- Helps to reduce wetness
- Controls odor
- Alcohol-free
- Light, fresh scent
- Made in USA



Item No.	Description	Quantity
3306PC5015	Roll-On Antiperspirant and Deodorant, 1.5 oz., Lightly scented	1 / EA, 96 / CA



PERSONAL CARE
Shampoo & Body Bath

- Gentle non-irritating formula
- Can be used on hair and body
- Convenient flip-top cap for easy dispensing
- pH balanced
- Pleasant and mild fragrance
- Made in USA



Item No.	Description	Quantity
3306PC8020	Shampoo and Body Bath, 2 oz. w/ Flip-top cap	1 / EA, 144 / CA



PILLOWS

Disposable Pillows

- Soft and comfortable spunbond polypropylene cover
- 100% garneted polyester fiber fill
- Non-allergenic, non-irritating, non-linting
- Single patient use
- Made in USA
- Compressed packaging saves storage space



Item No.	Description	Quantity
1233089EC100	Standard Disposable Pillow, White, 18" x 24"	1 / EA, 12 / CA
1233089DP10	Disposable Cough Pillow, White, 13" x 17", 8 oz.	1 / EA, 24 / CA
1233089DP20	Disposable Pillow, White, 16" x 22", 10 oz.	1 / EA, 12 / CA
1233089DP30	Disposable Pillow, White, 18" x 24", 12 oz.	1 / EA, 12 / CA
1233089DP40	Disposable Pillow, White, 18" x 24", 15 oz.	1 / EA, 12 / CA
1233089DP50	Disposable Pillow, White, 20" x 26", 17 oz.	1 / EA, 12 / CA
1233089DP55	Disposable Pillow, White, 18" x 24" 20 oz.	1 / EA, 12 / CA
1233089DP60	Disposable Pillow, White, 20" x 26", 20 oz.	1 / EA, 12 / CA
1233089DP65	Disposable Pillow, White, 18" x 24", 24 oz.	1 / EA, 12 / CA

Complementary Products

Item No.	Description
6358DPC001	Standard Disposable Pillowcase, White, 21" x 30"
6358DPC002	Premium Disposable Pillowcase, White, 21" x 30"
1233089LP30	Limited Reusable Pillow, Blue, 20" x 26", 18oz.



PILLOWS

Limited Reusable Pillows

- Cost-effective limited reuse pillow
- Moisture and stain resistant
- Multiple patient uses
- Made in USA



Item No.	Description	Quantity
1233089LP20	Limited Reusable Pillow, Blue, 18" x 24"	1 / EA, 12 / CA
1233089LP30	Limited Reusable Pillow, Blue, 20" x 26", 18 oz.	1 / EA, 12 / CA

Complementary Products

Item No.	Description
6358DPC001	Standard Disposable Pillowcase, White, 21" x 30"
6358DPC002	Premium Disposable Pillowcase, White, 21" x 30"



PILLOWS

Reusable Pillows

- Soft and comfortable Endurance™ vinyl ticking wipes clean*
- Model 1233089RP50 has easy-clean nylon ticking
- 100% garneted polyester fiber fill
- Fluid proof, flame retardant
- Anti-microbial, anti-fungal
- Non-allergenic, non-irritating, non-linting
- Made in USA
- Compressed packaging to save space

*Endurance™ is a trademark of Spec-Text, Inc.



Item No.	Description	Quantity
1233089RP10	Reusable Pillow, Blue, 13" x 17", 7 oz.	1 / EA, 24 / CA
1233089RP20	Reusable Pillow, Blue, 18" x 24", 17 oz.	1 / EA, 12 / CA
1233089RP30	Reusable Pillow, Blue, 20" x 26", 19 oz.	1 / EA, 12 / CA
1233089RP40	Reusable Pillow, Blue, 20" x 26", 24 oz.	1 / EA, 12 / CA
1233089RP40W	Reusable Pillow, White, 20" x 26", 24 oz.	1 / EA, 12 / CA
1233089RP50	Reusable Pillow, Yellow Nylon, 21" x 27", 22 oz.	1 / EA, 12 / CA
1233089RP60	Reusable Pillow, Yellow Nylon, 20" x 26", 18 oz.	1 / EA, 12 / CA
1233089RP65	Reusable Pillow, Yellow Nylon, 18" x 24", 16 oz.	1 / EA, 12 / CA

Complementary Products

Item No.	Description
6358DPC001	Standard Disposable Pillowcase, White, 21" x 30"
6358DPC002	Premium Disposable Pillowcase, White, 21" x 30".



Berman Airways

- Satin-smooth finish is safe and comfortable
- Size marked on bite block for easy identification
- Available in seven sizes, packaged in bulk or individually
- Manufactured in polypropylene



Item No.	Description	Quantity
0205BA1240	Berman Airway, 40mm, 12 per bag	1 / EA, 48 / CA
0205BA1250	Berman Airway, 50mm, 12 per bag	1 / EA, 48 / CA
0205BA1260	Berman Airway, 60mm, 12 per bag	1 / EA, 48 / CA
0205BA1270	Berman Airway, 70mm, 12 per bag	1 / EA, 48 / CA
0205BA1280	Berman Airway, 80mm, 12 per bag	1 / EA, 48 / CA
0205BA1290	Berman Airway, 90mm, 12 per bag	1 / EA, 48 / CA
0205BA12100	Berman Airway, 100mm, 12 per bag	1 / EA, 48 / CA
0205BA5040	Berman Airway, 40mm, Individually Bagged	1 / EA, 50 / CA
0205BA5050	Berman Airway, 50mm, Individually Bagged	1 / EA, 50 / CA
0205BA5060	Berman Airway, 60mm, Individually Bagged	1 / EA, 50 / CA
0205BA5070	Berman Airway, 70mm, Individually Bagged	1 / EA, 50 / CA
0205BA5080	Berman Airway, 80mm, Individually Bagged	1 / EA, 50 / CA
0205BA5090	Berman Airway, 90mm, Individually Bagged	1 / EA, 50 / CA
0205BA50100	Berman Airway, 100mm, Individually Bagged	1 / EA, 50 / CA



Mesh Briefs

- Soft, breathable mesh
- Now available in sizes Large to 4XL
- Color-coded stitching indicates size
- Available in bulk, 2 ea/bag and 5 ea/bag
- Made in USA



Item No.	Description	Quantity
4350MB4002	Mesh Brief, Large, 2/PK, Brown, 2/PK, 50 PK/CS	2 / BG, 100 / CA
4350MB4005	Mesh Brief, Large, 5/PK, Brown, 5/PK, 20 PK/CS	5 / BG, 100 / CA
4350MB4100	Mesh Brief, Large, Bulk, Brown, 100 EA/CS	1 / EA, 100 / CA
4350MB5002	Mesh Brief, XL, 2/PK, Green, 2/PK, 50 PK/CS	2 / BG, 100 / CA
4350MB5005	Mesh Brief, XL, 5/PK, Green, 5/PK, 20 PK/CS	5 / BG, 100 / CA
4350MB5100	Mesh Brief, XL, Bulk, Green, 100 EA/CS	1 / EA, 100 / CA
4350MB6002	Mesh Brief, XXL, 2/PK, Purple, 2/PK, 50 PK/CS	2 / BG, 100 / CA
4350MB6005	Mesh Brief, XXL, 5/PK, Purple, 5/PK, 20 PK/CS	5 / BG, 100 / CA
4350MB6100	Mesh Brief, XXL, Bulk, Purple, 100 EA/CS	1 / EA, 100 / CA
4350MB7002	Mesh Brief, 3XL, 2/PK, Lt. Blue, 2/PK, 50 PK/CS	2 / BG, 100 / CA
4350MB7005	Mesh Brief, 3XL, 5/PK, Lt. Blue, 5/PK, 20 PK/CS	5 / BG, 100 / CA
4350MB7100	Mesh Brief, 3XL, Bulk, Lt. Blue, 100 EA/CS	1 / EA, 100 / CA
4350MB8002	Mesh Brief, 4XL, 2/PK, White, 2/PK, 50 PK/CS	2 / BG, 100 / CA
4350MB8005	Mesh Brief, 4XL, 5/PK, White, 5/PK, 20 PK/CS	5 / BG, 100 / CA
4350MB8100	Mesh Brief, 4XL, Bulk, White, 100 EA/CS	1 / EA, 100 / CA



SLIPPERS Fall Management Slippers

- Bright colors visually alert staff to fall risk patients
- Sizes Small to Bariatric
- Double tread (top and bottom) to reduce fall risk if slipper rotates on the foot
- Large coverage, skid resistant, puff ink tread
- Soft terry cloth construction for optimal comfort
- Made in USA



Item No.	Description	Quantity
1296SLP7DR	Slipper, Red Fall Management, Small Double Tread, Smooth Side Out	2 / PR, 96 / CA
1296SLP7DY	Slipper, Yellow Fall Management, Small Double Tread, Smooth Side Out	2 / PR, 96 / CA
1296SLP7DG	Slipper, Green Fall Management, Small Double Tread, Smooth Side Out	2 / PR, 96 / CA
1296SLP8DR	Slipper, Red Fall Management, Medium Double Tread, Smooth Side Out	2 / PR, 96 / CA
1296SLP8DY	Slipper, Yellow Fall Management, Medium Double Tread, Smooth Side Out	2 / PR, 96 / CA
1296SLP8DG	Slipper, Green Fall Management, Medium Double Tread, Smooth Side Out	2 / PR, 96 / CA
1296SLP12DR	Slipper, Red Fall Management, Large Double Tread, Smooth Side Out	2 / PR, 96 / CA
1296SLP12DY	Slipper, Yellow Fall Management, Large Double Tread, Smooth Side Out	2 / PR, 96 / CA
1296SLP12DG	Slipper, Green Fall Management, Large Double Tread, Smooth Side Out	2 / PR, 96 / CA
1296SLP14DR	Slipper, Red Fall Management, XL Double Tread, Smooth Side Out	2 / PR, 96 / CA
1296SLP14DY	Slipper, Yellow Fall Management, XL Double Tread, Smooth Side Out	2 / PR, 96 / CA
1296SLP14DG	Slipper, Green Fall Management, XL Double Tread, Smooth Side Out	2 / PR, 96 / CA
1296SLP16DR	Slipper, Red Fall Management, XXL Double Tread, Smooth Side Out	2 / PR, 96 / CA
1296SLP16DY	Slipper, Yellow Fall Management, XXL Double Tread, Smooth Side Out	2 / PR, 96 / CA
1296SLP16DG	Slipper, Green Fall Management, XXL Double Tread, Smooth Side Out	2 / PR, 96 / CA
1296SLP18DR	Slipper, Red Fall Management, Bariatric Double Tread, Smooth Side Out	2 / PR, 96 / CA
1296SLP18DY	Slipper, Yellow Fall Management, Bariatric Double Tread, Smooth Side Out	2 / PR, 96 / CA
1296SLP18DG	Slipper, Green Fall Management, Bariatric Double Tread, Smooth Side Out	2 / PR, 96 / CA



SLIPPERS
Smooth Side Out Slippers

- Soft terry pile side of the slipper sock is on the inside for patient comfort and safety
- Large coverage, skid resistant, puff ink tread
- Available with double tread (top and bottom) to ensure grip regardless of how slipper is worn
- Color-coded by size
- Made in USA



Item No.	Description	Quantity
1296SLP1001	Slipper, Smooth Side Out, Single Tread, Toddler, White	2 / PR, 96 / CA
1296SLP1002	Slipper, Smooth Side Out, Single Tread, Small/Child, Yellow	2 / PR, 96 / CA
1296SLP1003	Slipper, Smooth Side Out, Single Tread, Medium, Blue	2 / PR, 96 / CA
1296SLP1004	Slipper, Smooth Side Out, Single Tread, Large, Beige	2 / PR, 96 / CA
1296SLP1005	Slipper, Smooth Side Out, Single Tread, XL, Gray	2 / PR, 96 / CA
1296SLP1006	Slipper, Smooth Side Out, Single Tread, XXL, Blue	2 / PR, 96 / CA
1296SLP2001	Slipper, Smooth Side Out, Double Tread, Toddler, White	2 / PR, 96 / CA
1296SLP2002	Slipper, Smooth Side Out, Double Tread, Small/ Child, Yellow	2 / PR, 96 / CA
1296SLP2003	Slipper, Smooth Side Out, Double Tread, Medium, Blue	2 / PR, 96 / CA
1296SLP2004	Slipper, Smooth Side Out, Double Tread, Large, Beige	2 / PR, 96 / CA
1296SLP2005	Slipper, Smooth Side Out, Double Tread, XL, Gray	2 / PR, 96 / CA
1296SLP2006	Slipper, Smooth Side Out, Double Tread, XXL, Blue	2 / PR, 96 / CA
1296SLP2007	Slipper, Smooth Side Out, Double Tread, Bariatric	2 / PR, 96 / CA
1296SLP4004	Slipper, Smooth Side Out, SingleTread, Large, Red	2 / PR, 96 / CA
1296SLP4005	Slipper, Smooth Side Out, Single Tread, Red, XL	2 / PR, 96 / CS



SLIPPERS
Textured Side Out Slippers

- Soft terry construction
- Skid-resistant, large-coverage treads
- Made in USA
- Color-coded by size
- Available in double tread (top and bottom) to ensure grip if sock rotates on the foot



Item No.	Description	Quantity
1296066699	Slipper, Terry, Non-slip Sole, Infant/Toddler, White	2 / PR, 96 / CA
1296066600	Slipper, Terry, Non-Skid Sole, Small/Child, Yellow	2 / PR, 96 / CA
1296066601	Slipper, Terry, Non-Skid Sole, Medium, Blue	2 / PR, 96 / CA
1296066602	Slipper, Terry, Non-Skid Sole, Large, Beige	2 / PR, 96 / CA
1296066603	Slipper, Terry, Non-Skid Sole, XL, Gray	2 / PR, 96 / CA
1296066604	Slipper, Terry, Non-Skid Sole, XXL, Lt. Blue	2 / PR, 96 / CA
1296066608	Slipper, Terry, Double Tread, Small/ Child, Yellow	2 / PR, 96 / CA
1296066609	Slipper, Terry, Double Tread, Medium, Blue	2 / PR, 96 / CA
1296066605	Slipper, Terry, Double Tread, Large, Beige	2 / PR, 96 / CA
1296066607	Slipper, Terry, Double Tread, XL, Gray	2 / PR, 96 / CA
1296066610	Slipper, Terry, Double Tread, XXL, Lt. Blue	2 / PR, 96 / CA
1296066611	Slipper, Terry, Double Tread, Bariatric, Beige	2 / PR, 96 / CA



TELEMETRY POUCHES Clear Polyethylene Telemetry Pouches

- Clear polyethylene pouch offers convenience, patient comfort, and protection for the transmitter
- 5¼" x 8⅞" pouch with clip helps protect transmitter from fluids
- Velcro® hook and loop closure easily allows transmitter wires to pass through top of pouch
- Pre-attached patient ID label



Item No.	Description	Quantity
35960MC500	Telemetry Pouch, 5¼" x 8⅞", Clear Polyethylene, Velcro® Hook & Loop Closure	1 / EA, 50 / CA



TELEMETRY POUCHES Foam Telemetry Pouches

- Foam pouch offers convenience, patient comfort, and protection for the transmitter
- Soft foam is comfortable for patient
- 4¾" x 6" fits most standard-size monitors
- Ideal for small digital "beeper"-size monitors
- Pre-attached patient ID label



Item No.	Description	Quantity
35960MC902	Telemetry Pouch, 4¾" x 6", Foam, Velcro® Hook & Loop, Waistbelt	1 / EA, 50 / CA



TELEMETRY POUCHES Non-Woven Telemetry Pouches with Neck and Waist Straps

- Non-woven pouch offers convenience, patient comfort, and protection for the transmitter
- Straps that tie at the neck and waist
- 6" x 9¼" pouch with non-woven back and polyethylene front
- Velcro® hook and loop closure
- Pre-attached patient ID label



Item No.	Description	Quantity
35960MC903	Telemetry Pouch, 6" x 9¼", Nonwoven Velcro® Hook & Loop, with Neck and Waist Straps	1 / EA, 50 / CA



TELEMETRY POUCHES Non-Woven Telemetry Pouches with Neck Straps

- Non-woven pouch offers convenience, patient comfort, and protection for the transmitter
- Straps that tie at the neck only
- 6" x 9¼" pouch with non-woven back and polyethylene front
- Velcro® hook and loop closure
- Pre-attached patient ID label



Item No.	Description	Quantity
35960MC901	Telemetry Pouch 6" x 9¼", Nonwoven, Velcro® Hook & Loop Closure w/ Neck Straps	1 / EA, 50 / CA



FACE MASKS
Procedure Masks

- Comfortable, cool and breathable
- Stretch earloops provide a nonbinding, personalized fit
- Convenient donning and quick fitting
- Fluid resistant to 80mm Hg
- Bacterial Filtration Efficiency > 99%
- Single use



Item No.	Description	Quantity
0205ELM001	Procedure Mask, Earloop, Yellow , 50 EA/BX, 10 BX/CS	1 / EA, 500 / CA



ISOLATION GOWNS
Cover Gowns

- Unrated for fluid protection; use when no fluid contact is expected
- Double poly spunbond, breathable, 23 GSM material
- Ties at neck and waist for proper fit
- Full back and elastic cuffs
- Good gown for visitors to isolation precaution rooms



Item No.	Description	Quantity
0205077750	Cover Gown, Unrated, Yellow Spunbond Universal Size, Full Back, Elastic Cuffs	1 / EA, 100 / CA
020577750XL	Cover Gown, Unrated, Yellow Spunbond XL Size, Full Back, Elastic Cuffs	1 / EA, 100 / CA



Personal Protective Gowns

- Use when heavy fluid contact is likely
- 0205077755 and 020577755XL meet ASTM F1671 for viral penetration
- Blue Polyethylene film provides high-level fluid protection
- Comfortable open-back design
- Thumbloop cuffs keep sleeves inside gloves for extra protection
- Universal and XL sizes
- Individually wrapped; 15 gowns per dispenser box
- Waist ties and apron neck with perforations for quick removal
- Glued seams



Item No.	Description	Quantity
0205077755	Protective Gown, Blue Polyethylene, Universal Open Back, Thumbloop Cuffs	1 / EA, 75 / CA
020577755XL	Protective Gown, Blue Polyethylene, XL, Open Back, Thumbloop Cuffs	1 / EA, 75 / CA



ISOLATION GOWNS Level 1 Isolation Gowns

- AAMI Level 1 Isolation Gown: Meets AAMI PB70 Level 1 testing
- Use when light fluid contact is expected
- Choose from tie or tape-tab neck closure
- Comfortable 22 GSM yellow SMS material
- Available in Universal and XL sizes
- Full back, ultrasonically welded seams and waist ties



Item No.	Description	Quantity
0205077751	Level 1 Isolation Gown, Tie Neck and Waist, Universal, Full Back, Elastic Wrist, 22 GSM Yellow SMS	1 / EA, 100 / CA
020577751XL	Level 1 Isolation Gown, Tie Neck and Waist, XL, Full Back, Elastic Wrist, 22 GSM Yellow SMS	1 / EA, 100 / CA



ISOLATION GOWNS Level 2 Isolation Gowns

- AAMI Level 2 isolation gown: Meets AAMI PB70 Level 2 testing
- Use when moderate fluid contact is possible
- Choose from tie or tape-tab neck closure
- Comfortable yet sturdy 25 GSM SMS material
- Available in Universal and XL sizes
- Ultrasonically welded seams
- Full back, waist ties, elastic cuffs



Item No.	Description	Quantity
0205077752	Level 2 Isolation Gown, Tie Neck and Waist, Universal, Full Back, Elastic Wrist, 25 GSM Yellow SMS	1 / EA, 100 / CA
020577752XL	Level 2 Isolation Gown, Tie Neck and Waist, XL, Full Back, Elastic Wrist, 25 GSM Yellow SMS	1 / EA, 100 / CA



ISOLATION GOWNS Level 3 Isolation Gowns

- AAMI Level 3 Isolation Gown
- Meets AAMI PB70 Level 3 testing
- Use when moderate to heavy fluid contact is possible
- Strong 36 GSM yellow SMS material
- Ultrasonically welded seams
- Available in Universal and XL sizes
- Full back, waist and neck ties, elastic cuffs



Item No.	Description	Quantity
0205077753	Level 3 Isolation Gown, Yellow SMS, Universal, Full Back, Elastic Wrist, Tie Neck and Waist	1 / EA, 100 / CA
020577753XL	Level 3 Isolation Gown, Yellow SMS, XL, Full Back, Elastic Wrist, Tie Neck and Waist	1 / EA, 100 / CA



DISPOSABLE COVERALLS Standard Disposable Coveralls

- Durable polypropylene material
- Knit cuffs
- Open collar and legs
- Full zip front
- Sizes Large (universal) to 3XL



Item No.	Description	Quantity
4363055803L	Coverall, Standard Polypropylene, Large, Knit Cuffs, Open Collar and Legs, Zipper Front	1 / EA, 24 / CA
4363055803X	Coverall, Standard Polypropylene, XL, Knit Cuffs, Open Collar and Legs, Zipper Front	1 / EA, 24 / CA
436355803XX	Coverall, Standard Polypropylene, 2XL, Knit Cuffs, Open Collar and Legs, Zipper Front	1 / EA, 24 / CA
436355803XXX	Coverall, Standard Polypropylene, 3XL, Knit Cuffs, Open Collar and Cuffs, Zipper Front	1 / EA, 24 / CA



DISPOSABLE COVERALLS Premium Disposable Coveralls

- Anti-static to <10kV for static discharge, static jump and static decay
- Made of 3-layer blue or white SMS material
- Sizes Large (universal) to 4XL
- Knit cuffs, full zip front, open collar and legs
- Meets OSHA Standard 29 CFR Part 1910.1030(d)(3)(xi)



Item No.	Description	Quantity
436305503L	Coverall, Premium, Anti-Static White SMS, Large, Knit Cuffs, Open Collar and Legs, Zip Front	1 / EA, 24 / CA
43635503BL	Coverall, Premium, Anti-Static Blue SMS, Large, Knit Cuffs, Open Collar and Legs, Zip Front	1 / EA, 24 / CA
436305503X	Coverall, Premium, Anti-Static White SMS, XL, Knit Cuffs, Open Collar and Legs, Zip Front	1 / EA, 24 / CA
43635503BX	Coverall, Premium, Anti-Static Blue SMS, XL, Knit Cuffs, Open Collar and Legs, Zip Front	1 / EA, 24 / CA
43635503XX	Coverall, Premium, Anti-Static White SMS, 2XL, Knit Cuffs, Open Collar and Legs, Zip Front	1 / EA, 24 / CA
43635503BXX	Coverall, Premium, Anti-Static Blue SMS, 2XL, Knit Cuffs, Open Collar and Legs, Zip Front	1 / EA, 24 / CA
43635503XXX	Coverall, Premium, Anti-Static White SMS, 3XL, Knit Cuffs, Open Collar and Legs, Zip Front	1 / EA, 24 / CA
43635503BXXX	Coverall, Premium, Anti-Static Blue SMS, 3XL, Knit Cuffs, Open Collar and Legs, Zip Front	1 / EA, 24 / CA
43635503XXXX	Coverall, Premium, Anti-Static White SMS, 4XL, Knit Cuffs, Open Collar and Legs, Zip Front	1 / EA, 24 / CA
4363503BXXXX	Coverall, Premium, Anti-Static Blue SMS, 4XL, Knit Cuffs, Open Collar and Legs, Zip Front	1 / EA, 24 / CA



DISPOSABLE LAB COATS Standard Disposable Lab Coats

- Made of sturdy polypropylene nonwoven material
- Open collar and knit cuffs
- 3 pockets, 5-snap front
- Sizes Small to 2XL
- White only



Item No.	Description	Quantity
436303304S	Lab Coat, Standard, White, Small, Snap Front, Knit Cuff	1 / EA, 25 / CA
436303304M	Lab Coat, Polypropylene, Medium, 3 Pockets, Knit Cuffs, Snap Front	1 / EA, 25 / CA
436303304L	Lab Coat, Polypropylene, Large, 3 Pockets, Knit Cuffs, Snap Front	1 / EA, 25 / CA
436303304X	Lab Coat, Polypropylene, XL, 3 Pockets, Knit Cuffs, Snap Front	1 / EA, 25 / CA
43633304XX	Lab Coat, Polypropylene, 2XL, 3 Pockets, Knit Cuffs, Snap Front	1 / EA, 25 / CA



DISPOSABLE LAB COATS Premium Disposable Lab Coats

- Antistatic to <10kv for static discharge, static jump and static decay
- Made of fluid-resistant SMS material
- Two color options: blue and white
- Meets OSHA Standard 29 CFR Part 1910.1030(d)(3)(xi)
- Open collar and knit cuffs
- 3 pockets, 5-snap front
- Sizes Small to 3XL



Item No.	Description	Quantity
436303504BS	Lab Coat, Anti-Static Blue SMS, Small, 3 Pockets, Knit Cuffs, Snap Front	1 / EA, 25 / CA
436303504BM	Lab Coat, Anti-Static Blue SMS, Medium, 3 Pockets, Knit Cuffs, Snap Front	1 / EA, 25 / CA
436303504BL	Lab Coat, Anti-Static Blue SMS, Large, 3 Pockets, Knit Cuffs, Snap Front	1 / EA, 25 / CA
436303504BX	Lab Coat, Anti-Static Blue SMS, XL, 3 Pockets, Knit Cuffs, Snap Front	1 / EA, 25 / CA
43633504BXX	Lab Coat, Anti-Static Blue SMS, 2XL, 3 Pockets, Knit Cuffs, Snap Front	1 / EA, 25 / CA
43633504BXXX	Lab Coat, Anti-Static Blue SMS, 3XL, 3 Pockets, Knit Cuffs, Snap Front	1 / EA, 25 / CA
436303504S	Lab Coat, Anti-Static White SMS, Small, 3 Pockets, Knit Cuffs, Snap Front	1 / EA, 25 / CA
436303504M	Lab Coat, Anti-Static White SMS, Medium, 3 Pockets, Knit Cuffs, Snap Front	1 / EA, 25 / CA
436303504L	Lab Coat, Anti-Static White SMS, Large, 3 Pockets, Knit Cuffs, Snap Front	1 / EA, 25 / CA
436303504X	Lab Coat, Anti-Static White SMS, XL, 3 Pockets, Knit Cuffs, Snap Front	1 / EA, 25 / CA
43633504XX	Lab Coat, Anti-Static White SMS, 2XL, 3 Pockets, Knit Cuffs, Snap Front	1 / EA, 25 / CA
43633504XXX	Lab Coat, Anti-Static White SMS, 3XL, 3 Pockets, Knit Cuffs, Snap Front	1 / EA, 25 / CA



Enzyme Gel Spray Pre-Cleaner

- Initiates enzymatic pre-cleaning
- Reduces heavy instrument transport weight compared to liquid-filled trays
- Free-rinsing and compatible with all washers
- Minimizes potential spill exposures due to splashing
- 32-oz. trigger bottle



Item No.	Description	Quantity
4067EGS100	Enzymatic Gel Spray Pre-Cleaner, 32 oz.	1 / EA, 6 / CA



Instrument Lubricants

- Neutral pH
- Contains anti-microbial ingredients
- Extends the life of the instruments and ensures proper operation of hinges and joints
- Easier cleaning: no rinsing or wiping
- Reduces rusting, discoloration, and corrosion
- Water soluble and steam permeable



Item No.	Description	Quantity
4067INL001	Instrument Lubricant, 1 Gallon	1 / EA, 4 / CA
4067INL015	Instrument Lubricant, 15 Gallon	1 / DR, 1 / EA



Enzyme Cleaner Pre-Soak

- Neutral pH, non-corrosive
- Three enzymes work together to remove soils quickly
- Protease: breaks down proteins
- Lipase: breaks down fat
- Amylase: breaks down starch and carbohydrates
- Disintegrates proteinaceous matter in hard-to-reach crevices



Item No.	Description	Quantity
4067ENZ001	Enzymatic Cleaner, 1 Gallon	1 / EA, 4 / CA
4067ENZ015	Enzymatic Cleaner, 15 Gallon	1 / EA, 1 / DR
4067MCPUN1G	Manual Pump for Enzymatic Cleaner	1 / EA
4067PEZ001	Premium Enzymatic Cleaner, 1 Gallon	1 / EA, 4 / CA
4067PEZ015	Premium Enzymatic Cleaner, 15 Gallon	1 / DR, 1 / EA



Neutral PH Instrument Cleaners

- Neutral pH
- Effective and safe way to maximize instrument cleanliness and service life
- Safe to use on brass, copper and aluminum
- Contains no phosphates
- Non-corrosive
- Rinse aid additive



Item No.	Description	Quantity
4067NPH001	Neutral PH Detergent, 1 Gallon	1 / EA, 4 / CA
4067NPH015	Neutral PH Detergent, 15 Gallon	1 / DR, 1 / EA



INSTRUMENTS
Thumb Forcep and Littauer Scissors

- High-quality, surgical-grade, stainless steel, single-use instruments
- Designed for cutting skin sutures
- One blade is hook-shaped to easily slip underneath suture to prevent it from slipping off the end of the blade
- Eliminates time and expense of reprocessing soiled instruments
- Meets ASTM and FDA requirements for surgical-grade steel



Item No.	Description	Quantity
6422FG0215	Littauer Scissor & Thumb Forcep, Sterile	1 / EA, 50 / CA
6422FG0213	Scissor, Littauer, Wire, Sterile, 4 1/2 "	1 / EA, 50 / CA
6422FG0214	Forcep, Thumb, Wire, Sterile, 4"	1 / EA, 50 / CA



Lubricating Jelly

- Intended for medical purposes to lubricate body orifices in facilitating the entry of diagnostic and therapeutic devices
- Provides instant lubricating action for patient comfort
- Non-alcohol and bacteriostatic
- Greaseless, water soluble and will not harm rubber appliances or medical instruments
- Available in 2-oz. and 4-oz. flip-top tubes or 3-gram and 5-gram foil packets



Item No.	Description	Quantity
4067063021	Lubricating Jelly, Sterile, 2 oz., Flip-top Tube	1 / EA, 72 / CA
4067063031	Lubricating Jelly, Sterile, 3 gram, Foil Packet	1 / EA, 1728 / CA
4067063041	Lubricating Jelly, Sterile, 4 oz., Flip-Top Tube	1 / EA, 72 / CA
4067063051	Lubricating Jelly, Sterile, 5 gram Foil Packet	1 / EA, 864 / CA



HAND-HELD SPHYGMOMANOMETERS Economy Aneroid Sphygmomanometers

- Standard, chrome-plated deflation valve
- Latex inflation bladder, bulb and tubing
- 10-year calibration warranty on manometer
- Common features to Economy Model:
 - Black enamel, 300mmHg, no-pin-stop, spring-loaded manometer
 - Measures blood pressure from 20mmHg to 300mmHg, certified accurate to +/-3mmHg



Item No.	Description	Quantity
3596994100	Sphygmomanometer, Adult Economy	1 / EA, 1 / BX



HAND-HELD SPHYGMOMANOMETERS
Latex-Free Economy Aneroid Sphygmomanometers

- Standard, chrome-plated deflation valve
- Dip-molded inflation bladder withstands heavy use
- 10-year calibration warranty on manometer
- Common features to Economy Model:
 - Black enamel, 300mmHg, no-pin-stop, spring-loaded manometer
 - Measures blood pressure from 20mmHg to 300mmHg, certified accurate to +/-3mmHg



Item No.	Description	Quantity
3596994105	Sphygmomanometer, Adult Economy, Latex-Free	1 / EA, 1 / BX



HAND-HELD SPHYGMOMANOMETERS
Deluxe Aneroid Sphygmomanometers

- Latex inflation bulb fitted with filter-screen to reduce dust build-up and ensure precise deflation control and reliable service
- Dip-molded, latex inflation bladder withstands heavy use
- Lifetime calibration warranty on manometer
- Common features to Deluxe Model:
 - Black enamel, 300mmHg, no-pin-stop, spring-loaded manometer
 - Measures blood pressure from 20mmHg to 300mmHg, certified accurate to +/-3mmHg



Item No.	Description	Quantity
3596994200	Sphygmomanometer Adult Deluxe	1 / EA, 1 / BX
3596994210	Sphygmomanometer Large Adult Deluxe	1 / EA, 1 / BX
3596994230	Sphygmomanometer Child Deluxe	1 / EA, 1 / BX



HAND-HELD SPHYGMOMANOMETERS
Latex-Free Deluxe Aneroid Sphygmomanometers

- Inflation bulb fitted with filter-screen to reduce dust build-up and ensure precise deflation control and reliable service
- Dip-molded inflation bladder withstands heavy use
- Lifetime calibration warranty on manometer
- Common features to Deluxe Model:
 - Black enamel, 300mmHg, no-pin-stop, spring-loaded manometer
 - Measures blood pressure from 20mmHg to 300mmHg, certified accurate to +/-3mmHg



Item No.	Description	Quantity
3596994205	Sphygmomanometer, Adult, Deluxe, Latex-Free	1 / EA, 1 / BX
3596994215	Sphygmomanometer, Large, Adult, Deluxe, Latex-Free	1 / EA, 1 / BX
3596994235	Sphygmomanometer, Child, Deluxe, Latex-Free	1 / EA, 1 / BX



STETHOSCOPES
Disposable Chrome and Plastic Stethoscopes

- Lightweight plastic, diaphragm-only blood pressure stethoscope with chrome-plated binaurals
- Recommended for single patient use
- PVC tubing in yellow or black
- Overall length: 30"; diaphragm diameter: 1¾"
- Weight: 2½ oz.
- Individually packaged in a sealed plastic bag



Item No.	Description	Quantity
3596994700	Stethoscope Disposable Chrome Black	1 / EA, 10 / BX
3596994730	Stethoscope Disposable Chrome Yellow	1 / EA, 10 / BX



STETHOSCOPES
Disposable Plastic Stethoscopes

- Lightweight plastic, diaphragm-only blood pressure stethoscope
- Recommended for single patient use
- PVC tubing in yellow or cobalt blue
- Overall length: 30"; diaphragm diameter: 1¾"
- Weight: 1½ oz.
- Individually packaged in a sealed plastic bag



contact your local Owens & Minor representative for details

Item No.	Description	Quantity
3596994910	Stethoscope Disposable Plastic Blue	1 / EA, 10 / BX
3596994930	Stethoscope Disposable Plastic Yellow	1 / EA, 10 / BX



STETHOSCOPES
Neonatal Dual-Head Stethoscopes

- Chrome-plated brass combination chest piece with 25mm (.975") diaphragm and bell
- Slip-on, non-chill bell ring
- Adjustable chrome-plated brass binaural with PVC ear tips
- Flexible black PVC tubing
- Overall length: 30"
- Weight: 4 oz.



Item No.	Description	Quantity
3596995100	Stethoscope Neonatal Dual-Head Black	1 / EA, 1 / BX



STETHOSCOPES Infant Dual-Head Stethoscopes

- Chrome plated zinc combination chest-piece with 25mm diaphragm and bell
- Slip-on, gray, non-chill bell ring for patient comfort
- Adjustable chrome-plated brass binaural with PVC ear tips
- Overall length 32"
- Weight: 4 oz.



Item No.	Description	Quantity
3596995250	Stethoscope Infant Dual-Head Pink	1 / EA, 30 / CA



STETHOSCOPES Pediatric Dual-Head Stethoscopes

- Lightweight aluminum combination chest piece with 1¼" non-chill bell with non-chill rim
- Adjustable chrome-plated brass binaural with PVC ear tips
- Flexible black PVC tubing
- Overall length: 30"
- Weight: 4 oz.



Item No.	Description	Quantity
3596995000	Stethoscope Pediatric Dual Head Black	1 / EA, 1 / BX
3596995020	Stethoscope Pediatric Dual Head Grey	1 / EA, 30 / CA



STETHOSCOPES Lightweight Single-Head Stethoscopes

- Lightweight aluminum, diaphragm-only stethoscope with nurse-style chest piece
- Adjustable chrome-plated brass binaural with PVC ear tips
- Flexible PVC tubing in black, red, blue, green, gray, gold, yellow, light blue, teal or white
- Overall length: 30"; diaphragm diameter: 1¼"
- Weight: 4 oz.



Item No.	Description	Quantity
3596994600	Stethoscope Single Head Black	1 / EA, 1 / BX
3596994605	Stethoscope Single Head Red	1 / EA, 1 / BX
3596994610	Stethoscope Single Head Blue	1 / EA, 1 / BX
3596994615	Stethoscope Single Head Green	1 / EA, 1 / BX
3596994620	Stethoscope Single Head Gray	1 / EA, 1 / BX
3596994625	Stethoscope Single Head Gold	1 / EA, 1 / BX
3596994630	Stethoscope Single Head Yellow	1 / EA, 1 / BX
3596994635	Stethoscope Single Head Teal	1 / EA, 30 / CA
3596994640	Stethoscope Single Head Lt Blue	1 / EA, 1 / BX
3596994645	Stethoscope Single Head White	1 / EA, 1 / BX



STETHOSCOPES
Lightweight Dual-Head Stethoscopes

- Lightweight aluminum combination diaphragm/bell chest piece for wide frequency response
- Non-chill 1 3/4" diaphragm and 3/4" bell
- Adjustable, chrome-plated brass binaural with PVC ear tips
- Flexible PVC tubing in black, red, blue, green, gray, teal, neon pink or yellow
- Overall length: 30"; diaphragm diameter 1 3/4"
- Weight: 4 oz.



Item No.	Description	Quantity
3596994300	Stethoscope Dual Head Black	1 / EA, 1 / BX
3596994305	Stethoscope Dual Head Red	1 / EA, 1 / BX
3596994310	Stethoscope Dual Head Blue	1 / EA, 1 / BX
3596994315	Stethoscope Dual Head Green	1 / EA, 1 / BX
3596994320	Stethoscope Dual Head Gray	1 / EA, 1 / BX
3596994330	Stethoscope Dual Head Yellow	1 / EA, 1 / BX
3596994335	Stethoscope Dual Head Teal	1 / EA, 30 / CA
3596994350	Stethoscope Dual Head Pink	1 / EA, 1 / BX

contact your local Owens & Minor representative for details



STETHOSCOPES
Lightweight Single-Head Bowles Stethoscopes

- Chrome-plated zinc, diaphragm-only chest piece with traditional, raised-stem design
- Adjustable, chrome-plated brass binaural with PVC ear tips
- Flexible, black PVC tubing
- Overall length: 30"; diaphragm diameter: 1.625"
- Weight: 4 oz.



Item No.	Description	Quantity
3596994800	Stethoscope Bowles Black	1 / EA, 1 / BX



STETHOSCOPES
Professional Adult and Pediatric Stethoscopes

- Full-range stethoscope with unmatched comfort and superior acoustics
- Comfort sealing eartips and lightweight design
- Superior acoustics allow for accurate auscultation of cardiac, pulmonary and vascular sounds
- Diaphragm is excellent for all sounds and murmurs. Best for higher frequencies and lung sounds. Ideal in identifying splitting of sounds, systolic clicks and ejection sounds
- Bell provides excellent detection and clarity for low and medium frequency heart sounds and murmurs such as a faint gallop or diastolic rumble
- Overall length: 28"; (Adult) 4.5 cm; bell diameter 3.0 cm (Pediatric) 3.5 cm; bell diameter 2.5 cm



Item No.	Description	Quantity
3596MC79135	Stethoscope Adult Professional Dual- Head Black	1 / EA
3596MC79137	Stethoscope Adult Professional Dual- Head Navy	1 / EA
3596MC79139	Stethoscope Adult Professional Dual-Head Burgundy	1 / EA
3596MC79145	Stethoscope Ped Professional Dual-Head Black	1 / EA



STETHOSCOPES Dual-Head Cardiology Stethoscopes

- Full-range cardiology/pulmonology stethoscope
- Can identify extremely low frequencies such as a faint diastolic gallop or rumble, as well as the higher frequencies of the faint early blowing aortic diastolic murmur
- Diaphragm is excellent for all sounds and murmurs. Best for higher frequencies such as faint diastolic murmur of aortic regurgitation. Ideal in identifying splitting of sounds, systolic clicks and ejection sounds
- Diaphragm is designed to accentuate sounds in the 200 to 500 Hz range—the range where many significant heart sounds occur
- Bell provides excellent detection and clarity for low and medium frequency heart sounds and murmurs (such as a faint gallop or diastolic rumble)
- Overall length: 25"; diaphragm diameter 4.85 cm; bell diameter 3.25 cm



Item No.	Description	Quantity
3596MC79125	Stethoscope Cardiology Elite Dual-Head Black	1 / EA



STETHOSCOPES Sprague-Rappaport Stethoscopes with Accessory Kit

- Threaded chest piece drum of chrome-plated zinc with interchangeable fittings and other accessories for complete diagnostic capability
- Double-tube configuration for better acoustic separation
- Valve mechanism eliminates acoustic leakage
- Adjustable chrome-plated brass binaural with PVC ear tips
- Accessory kit includes adult and pediatric transparent diaphragms; three bells (1¼", 1", ¾"); one pair mushroom ear tips; one pair silicon ear tips; and a chest piece key for easy diaphragm replacement
- Flexible PVC tubing available in black or neon yellow



Item No.	Description	Quantity
3596994510	Stethoscope Sprague-Rappaport 22" Blk	1 / EA, 1 / BX
3596994535	Stethoscope Sprague-Rappaport 22" NYL	1 / EA, 1 / BX



DIGITAL THERMOMETERS
Oral and Rectal Digital Thermometers

- Measures temperature in both Fahrenheit and Celsius: 89.6°–109.4°F, 32.0°–43.9°C (dual scale)
- Temperature reading can take approximately:
 - Oral: 40–90 seconds
 - Rectal: 35–60 seconds
 - Axillary: 60–180 seconds
- Automatic shut-off
- Liquid crystal display (LCD)
- Performs automatic self-test to ensure accuracy
- Microtechnology ensures precision measurement within $\pm 0.2^{\circ}\text{F}$ and $\pm 0.1^{\circ}\text{C}$



Item No.	Description	Quantity
3596916814A	Thermometer, Digital, Oral, Dual Scale, Blue Cap	1 / EA, 25 / CA
3596916815A	Thermometer, Digital, Rectal, Dual Scale, Red Cap	1 / EA, 25 / CA
3596916820A	Thermometer, Digital, Oral, Rectal, and Axillary Use, Dual Scale, Neutral Cap	1 / EA, 25 / CA



DIGITAL THERMOMETERS
Oral Thermometer Kits with Covers

- Microtechnology ensures precision measurement within $\pm 0.2^{\circ}\text{F}$ and $\pm 0.1^{\circ}\text{C}$
- Measures temperatures from 89.6–109.4°F and 32.0–43.9°C
- Temperature reading for oral thermometer takes 40–90 seconds
- Liquid crystal display (LCD)
- "On-off" switch and automatic switch-off
- Performs automatic self-test to ensure accuracy
- Sterile thermometer sheath for single patient use



Item No.	Description	Quantity
3596916812A	Thermometer, Digital, Oral, with 20 Thermometer Sheaths, Dual Scale, Blue Cap	1 / EA, 20 / CA
3596206335	Thermometer Sheath for Digital Thermometer	1 / EA, 5000 / CA



EKG Paper

- Chemical Thermal coated paper manufactured to original equipment manufacturer's specifications
- Consistent trace quality
- Fully warranted materials and workmanship
- Manufactured in the USA
- 5-year warranty



Item No.	Description	Quantity
4067M0007983	EKG Chart Paper, 3-Channel, Red Grid, 216mm x 280mm, Z-Fold	1 / EA, 10 / CA
4067M0997984	EKG Chart Paper, 3-Channel, Red Grid, 216mm x 280mm, Z-Fold	200 / PK, 2000 / CA
4067M0997989	EKG Chart Paper, 3-Channel, Red Grid, 216mm x 280mm, Z-Fold	1 / EA, 10 / CA
4067M4305AAO	EKG Fetal Chart Paper, Red Grid, Z-Fold	1 / EA, 40 / CA
4067M4305CAO	EKG Fetal Chart Paper, Red Grid, 152mm x 90mm, With Adhesive Labels, Z-Fold	1 / PK, 40 / CA
4067MC008038	EKG Chart Paper, Single Channel-Blank, 50mm x 100'	1 / RL, 100 / CA
4067MC1100HD	Viedo Print Paper, 110mm x 20mm	1 / RL, 5 / BX
4067MC18416	EKG Chart Paper, Red Grid, 3-Channel, Full Grid Stress Paper, 216mm x 280mm, Z-Fold	1 / PK, 20 / CA
4067MC2485A	EKG Chart Paper, Red Grid, 3-Channel, 216mm x 280mm, Archive Grade (25 years)	1 / EA, 8 / CA
4067MC40477	EKG Chart Paper, Single Channel, No Grid, (Machine Prints Grid), 50mm x 100'	1 / RL, 100 / CA
4067MC9420	EKG Chart Paper, Red Grid, Stress Paper, 3-Channel, 216mm x 280mm	1 / PK, 8 / CA
4067MC9424	EKG Chart Paper, Red Grid, 3-Channel With Header, 216mm x 280mm	1 / PK, 8 / CA

Item No.	Description	Quantity
4067MCM1707A	EKG Chart Paper, Red Grid, 3-Channel, 216mm x 280mm, Z-Fold	1 / PK, 8 / CA
4067MCM1910A	EKG Fetal Chart Paper, Orange Grid, 150mm x 100mm, Z-Fold	1 / PK, 40 / CA
4067MCMF2428	EKG Mounting Form For Single Channel Chart Paper, 50mm x 120mm, Single Sheet	100 / PK, 500 / CA
4067MCPCLP12	EKG Chart Paper, Single Channel, Red Grid, 106mm x 100'	1 / EA, 18 / CA
4067MSL60196	EKG Chart Paper, Single Channel, 50mm x 120mm, Z-Fold	1 / EA, 100 / CA



BLANKET WARMERS
Counter-Top Blanket Warmers

- Single-compartment, stainless steel blanket warmer with digital temperature control to ensure fast and easy operation
- Glass or solid self-closing, energy efficient door
- 1 fully adjustable perforated solid stainless steel shelf
- ETL approved and meets AORN and ECRI guidelines for patient safety
- Temperature range 90°–140°F; AORN and ECRI recommended temperature 120°F
- Power Supply: 120V, 50/60 cycles, and 5.8 amps
- Dimensions:
 - Overall Size: 30"W x 24"D x 24"H
 - Compartment Size: 28"W x 17½"D x 16"H (4.5 cubic feet)



Item No.	Description	Quantity
4347BW3024G	Counter-Top Blanket Warmer Glass Door	1 / EA
4347BW3024S	Counter-Top Blanket Warmer Solid Door	1 / EA



BLANKET WARMERS
Under-Counter Blanket Warmers

- Single-compartment stainless steel blanket warmer with digital temperature control to ensure fast and easy operation
- Glass or solid self-closing, energy efficient door
- 2 fully adjustable perforated solid stainless steel shelves
- ETL approved and meets AORN and ECRI guidelines for patient safety
- Temperature range 90°–140°F; AORN and ECRI recommended temperature 120°F
- Power Supply: 120V, 50/60 cycles, and 7.3 amps
- Dimensions:
 - Overall Size: 30"W x 24"D x 36"H
 - Compartment Size: 28"W x 17½"D x 30"H (8.45 cubic feet)



Item No.	Description	Quantity
4347BW3036G	Under-Counter Blanket Warmer Glass Door	1 / EA
4347BW3036S	Under-Counter Blanket Warmer Solid Door	1 / EA

Complementary Products

Item Number	Description
4347BW3504	Mobile Dolly, Stainless Steel, 4"H



BLANKET WARMERS Single-Door Blanket Warmers

- Single-compartment stainless steel blanket warmer with digital temperature control to ensure fast and easy operation
- Glass or solid self-closing, energy efficient door
- 3 fully adjustable perforated solid stainless steel shelves
- ETL approved and meets AORN and ECRI guidelines for patient safety
- Temperature range 90°–140° F; AORN and ECRI recommended temperature 120°F
- Power Supply: 120V, 50/60 cycles, and 7.3 amps
- Dimensions:
 - Overall Size: 30"W x 24"D x 72"H
 - Compartment Size: 28"W x 17½"D x 65"H (18.43 cubic feet)



Item No.	Description	Quantity
4347BW3072G	Single-Door Blanket Warmer Glass Door	1 / EA
4347BW3072S	Single-Door Blanket Warmer Solid Door	1 / EA
4347BW3504	Mobile Dolly Stainless Steel 4"H	1 / EA

Complementary Products

Item Number	Description
4347BW3504	Mobile Dolly, Stainless Steel, 4"H



BLANKET WARMERS Double-Door Blanket Warmers

- Dual-compartment stainless steel blanket warmer with independent digital temperature controls to ensure fast and easy operation
- Glass or solid self-closing, energy efficient door
- Upper and lower compartment each have 2 adjustable perforated solid stainless steel shelves
- ETL approved and meets AORN and ECRI guidelines for patient safety
- Temperature range 90°–140°F; AORN and ECRI recommended temperature 120°F
- Power Supply: 120V, 50/60 cycles, and 7.3 amps
- Dimensions:
 - Overall Size: 30"W x 24"D x 72"H
 - Compartment Size: 28"W x 17½"D x 30"H upper 35"H lower (18.43 total cubic feet)



Item No.	Description	Quantity
4347BW3072DG	Double-Door Blanket Warmer Glass Door	1 / EA
4347BW3072DS	Double-Door Blanket Warmer Solid Door	1 / EA

Complementary Products

Item Number	Description
4347BW3504	Mobile Dolly, Stainless Steel, 4"H



IRRIGATION TRAYS
Irrigation Trays with Bulb Syringe

- Easy-to-use syringes for irrigation or enteral feeding
- Soft pliable bulb w/ large finger flanges for one-handed operation and optimum control
- Easy-to-read graduations, 60cc capacity
- Top-quality components
- Complete offering of irrigation trays and components available to fit a variety of procedural needs



Item No.	Description	Quantity
0205ITRAY120	Irrigation Tray w/60cc Bulb Syringe w/ Tip Protector 1000cc Basin, 500cc Grad Container, & Patient ID Label	1 / EA, 20 / CA
0205ITRAY130	Irrigation Tray w/60cc Bulb Syringe w/ Tip Protector 1200cc Basin, 500cc Grad Container, Alcohol Prep Pad, Underpad, & Patient ID Label	1 / EA, 20 / CA
0205ITRAY140	CSR Wrapped Irrigation Tray w/60cc Bulb Syringe w/Tip Protector 1000cc Contoured Basin, 500cc Grad Container, Patient ID Label, & Waterproof CSR Wrap	1 / EA, 20 / CA

Complementary Products

Item Number	Description
0205BULB011P	Bulb Irrigation Syringe, 60ml, Latex-Free



IRRIGATION TRAYS
Irrigation Trays with Piston Syringe

- Easy-to-use syringes for irrigation or enteral feeding
- Ergonomic design with thumb control ring and large finger flanges for one-handed operation and optimum control
- Smooth-motion plunger with three-tiered, leak-proof seal
- Easy-to-read graduations, 60cc capacity
- Top-quality components
- Complete offering of irrigation trays and components available to fit a variety of procedural needs



Item No.	Description	Quantity
0205PTRAY126	Irrigation Tray w/60cc Piston Syringe w/ Tip Protector 1000cc Basin, 500cc Grad Container, & Patient ID Label	1 / EA, 20 / CA
0205PTRAY136	Irrigation Tray w/60cc Piston Syringe w/ Tip Protector 1200cc Basin, 500cc Grad Container, Alcohol Prep Pad, Underpad, & Patient ID Label	1 / EA, 20 / CA
0205PTRAY146	CSR Wrapped Irrigation Tray w/60cc Piston Syringe w/Tip Protector 1000cc Contoured Basin, 500cc Grad Container, Patient ID Label, & Waterproof CSR Wrap	1 / EA, 20 / CA

Complementary Products

Item Number	Description
0205PIST015P	Piston Irrigation Syringe, 60ml, Latex-Free



Incision and Drainage Trays

- Stackable tray, wrapped and bagged
- Most commonly used components
- Packaged in typical order of use



Central Line Dressing Trays

- Sterile trays protect antibiotic contents
- Stackable containers for safe, easy storage
- Convenient easy-to-open packaging
- Packaged in typical order of use



Incision & Drainage Trays

4366INC1001B

Waste Bag, 12" x 18", White w/Tie	1
Fenestrated Drape, 13" x 18"	1
Towel, 13" x 19", 2-ply, Tissue/Poly	1
Gauze, 3" x 3", 12-ply	4
PVP Prep Pad	1
Alcohol Prep Pad, Medium	1
Scalpel, #11	1
Forceps, 5", Wireform	1
Scissors, Littauer, 4½", Wireform	1
Hemostat, Kelly, Serrated, 5½"	1
Quantity/Case	20

Central Line Dressing Trays	4366CLC1001	4366CLC1002	4366CLC1004	4366CLC1007	6422CLC1014	6422CLC1015	4366CLC1017	4366CLC1023	4366CLC2001	4366CLC3001	6422DC1002	6422CLC5004	4366CLC1030
Packaging													
Soft Tray (Blister) w/ Tyvek® Lid	1	1	1	1	1	1	1	1	1	1		1	1
Hard Tray w/ Tyvek® Lid											1		
CSR Wrap 24" x 24"											1		
Prep and Dressing													
Gauze, 2" x 2", 8 Ply	3	3		1			2	3					2
Gauze, 4" x 4", 8 Ply	2		2					2		2		2	
Gauze, 2" x 2", 12 Ply											1		
Gauze 4" x 4", 12 Ply											1		
Gauze, Split, 2" x 2", Non-Woven		1											
Tegaderm™, 3½" x 4½" (#1655)				1					1				
Tegaderm™, 4" x 4¾" (#1616)	1	1	1		1			1			1	1	
Tegaderm™, 4" x 4¾" w/Frame (#16126)													1
Tegaderm™, CHG 3½" x 4" (#1657)							1						
Opsite® (#4016)									1				
MediChoice® Transparent Dressing 4" x 4½"						1							
Biopatch® Protective Disk w/CHG					1			1					1
Dressing, Non-Adherent , 2" x 2"											1		
Dressing, Non-Adherent , Split, 2" x 2"											1		
Towels/Drapes													
Drape, 13" x 19", 2 Ply, T/P, Flat	1	1	1				1	1	1	1		1	
Drape, 13" x 19½", 3 Ply, T/P/T				1									
Drape, 17" x 19"						1							
Antiseptics													
Alcohol Swabstick, 1/Pk											1		1
Alcohol Swabsticks, 3/Pk	1	1	1					1	1	1		1	
Skin Protectant Pad				1	1		1						
Skin Protectant Swab, 1/Pk		1											1
PVP Ointment		1											
PVP Swabsticks, 3/Pk	1	1						1					
ChloroPrep®, 3mL One Step w/ Insert	1		1	1			1	1	1	1	1		1
Chlorascrub™ Swabstick 1.6mL w/Insert					2							1	
Chlorascrub™ Swabstick 5.1mL w/Insert						1							
Miscellaneous													
StatLock®					1								
Transpore™ Tape, ¾" x 18"						1							
Transpore™ Tape, 1" x 18"	1	1						1	1		1		
Tape Strip, ½" x 4", 3/Cd			1				1					1	1
Cotton-Tip Applicator, Wood, 6"	1							1					
Dressing Change Label	1	1	1		2	1		1	1	1	1	1	
Gloves, Walleled, LF/PF, Medium Pair	1	1	1	1			1	1	1	1	1	1	1
Gloves, Nitrile, Walleled, Medium						2							
Gloves, Next Gen Vinyl, Cuffed, Pair					1								
Mask, Pleated Earloops, Blue	1	1	1				1	1	1	1			1
Mask, Ear Loop, Yellow						2					1	1	
Tape Measure							1						
Quantity/Case	30	30	30	30	30	30	30	30	30	30	20	30	30



IV START KITS IV Start Kits

- Stackable containers for safe, easy storage
- Convenient easy-to-open packaging
- Packaged in typical order of use



IV Start Kits	43661VS1001	43661VS1002	43661VS1006	43661VS1011	43661VS1012	43661VS1013	64221VS1052	64221VS1052B
Tray								
Soft Tray (Blister) w/ Tyvek® Lid	1	1	1	1	1	1	1	1
Dressings								
Gauze, 2" x 2", 8-Ply	2	2	2	2	2	2		
Gauze, 2" x 2", 4-Ply, Non-Woven							2	2
Tegaderm™ PF w/Border, 2½" x 2¾" (#1614)								1
Tegaderm™ PF, 2½" x 2¾" (#1624)	1				1	1		
OpSite® IV3000 Dressing (#4017D)		1	1					
Transparent Dressing, 2½" x 2¾"							1	
Towels/Drapes								
Wrap, 10" x 13", T/P						1		
Miscellaneous								
Transpore™ Tape, ¾" x 18"	1	1	1	1	1	1		1
Transparent Tape, ¾" x 18"							1	
Dressing Change Label	1	1	1	1	1	1	1	1
Tourniquet LF ¾" x 18" Rolled	1	1	1	1	1	1		
Tourniquet LF, 1" x 18"							1	
Tourniquet LF, 1" x 18", Orange								1
Adhesive Bandage				1				
Gloves, Exam, Latex-Free, Medium					1			
Antiseptics								
Alcohol Prep Pads	1	1		2	1	2	1	1
PVP Pad				1		1		
Tincture of Iodine, 2% (BD Persist)	1	1			1			
Quantity/Case	50	50	50	50	50	50	100	100



IV START KITS IV Start Kits with Chlorascrub

- Stackable containers for safe, easy storage
- Convenient easy-to-open packaging
- Packaged in typical order of use
- Includes Chlorascrub™ antiseptic



IV Start Kits with Chlorascrub

	43661VS1051	64221VS1001	64221VS8003	64221VS1053	64221VS1056	64221VS1059	64221VS1061A	43661VS1062	43661VS1066
Tray									
Soft Tray (Blister) w/ Tyvek® Lid	1	1	1	1	1	1	1	1	1
Dressings									
Gauze, 2" x 2", 8-Ply	2		2					2	
Gauze, 2" x 2", 4-Ply, Non-Woven		2		2	2	2	2		2
Tegaderm™ PF, 2¾" x 2¾" (#1624)	1		1			1			1
Tegaderm™ PF, w/ Border 2¾" x 2¾" (#1614)					1		1		
MediChoice® Transparent Dressing, 2¾" x 2¾"		1							
Transparent Dressing, 2¾" x 2¾"				1					
Towels/Drapes									
Wrap, 10" x 13", T/P								1	
Miscellaneous									
Transpore™ Tape, ¾" x 18"	1	1	1		1				1
Transparent Tape, ¾" x 18"				1			1		
Durapore™ Tape, 1" x 18"								1	
Dressing Change Label		1	1	1	1	1	1	1	1
Tourniquet LF, ¾" x 18", Rolled	1							1	1
Tourniquet LF, 1" x 18"		1	1	1	1	1	1		
Statlock® Securement Device						1			
Pre-Filled Syringe, Saline 5mL						1			
Antiseptics									
Chlorascrub™ Pad, 1.0 mL w/ Insert	1	1	1	1	1	1		1	1
Chlorascrub™ Short Swabstick, 1.6 mL w/ Insert							1		
Quantity/Case	50	50	50	100	50	50	50	50	50

Choose MediChoice®
Transparent Dressings
to protect I.V. sites and
surgical incisions.

Our dressings conform to hard-to-dress areas, allow constant observation, and can be left in place for several days.

Refer to the full MediChoice Catalog for additional item information.



IV START KITS IV Start Kits with ChloroPrep®

- Stackable containers for safe, easy storage
- Convenient easy-to-open packaging
- Packaged in typical order of use
- Includes ChloroPrep® antiseptic





Suture Trays

- Stackable tray, wrapped and bagged
- Most commonly used components
- Packaged in typical order of use
- Precision floor grade or economy instruments



Suture Trays

4366MSUT1002

6422MSUT1003

6422MSUT1004

	4366MSUT1002	6422MSUT1003	6422MSUT1004
Packaging			
Hard Tray, 1 Compartment	1		
Hard Tray, 2 Compartment			1
Hard Tray w/ Tyvek® Lid		1	
Wrap, CSR, 24" x 24"	1		
Prep and Dressing			
Gauze, 2" x 2", 8-ply	8		
Gauze, 4" x 4", 8-ply	5	10	
Gauze, 4" x 4", 12-ply			15
Instruments			
Webster, Needle Holder, Satin, 5"	1		1
Webster, Needle Holder, Smooth, 5"		1	
Scissors, Iris, Straight, 4½"	1	1	1
Forcep, Adson, Serrated, 4¾"		1	1
Forcep, Adson, Tissue, 4¾", 1 x 2 Teeth	1		
Forcep, Kelly, Straight 5½"			1
Antiseptics			
Lidocaine 2% 10mL			1
Neosporin®			1
Towels and Drapes			
Drape, Fenestrated, 18" x 26"	1	1	1
Towel, 17" x 19" Poly 2/Pk		1	
Towel, Kaycel®, Blue, 19" x 27"	1		
Towel, Kaycel®, White, 19" x 27"			2
Miscellaneous			
Safety Needle 25 GA x 5/8"			1
Safety Needle BD 18 GA x 1½"			1
Safety Needle BD 25 GA x 1½"			1
Syringe 12cc L/L			1
Medicine Cup, 2 oz.	2	1	
Ampule Holder 5 mL			1
Quantity/Case	20	20	20



Skin Staple Remover Trays

- Sharp instrument every time
- Easy, stable finger-control
- Fine point for easy access to staple
- Convenient easy-to-open packaging
- Packaged in typical order of use



Custom minor procedure kits & trays

You can pick and choose your own components including your preferred packaging, prep and dressing, towels and drapes, antiseptics, instruments and more.

- Custom kits built to meet your specific procedural requirements
- Choose from a wide range of brand-name components, as well as MediChoice® brand
- CHG Chloraprep® and Chlorascrub™ antiseptics, plus traditional PVP and alcohol
- Components can be arranged in order of use

Skin Staple Remover Trays	4366SSR1001	4366SSR1002	6422SSR1003
Packaging			
Hard Tray w/ Tyvek® Lid	1	1	1
Prep and Dressing			
Gauze, 3" x 3", 12-ply	1	1	1
Instruments			
Skin Staple Remover, Tweezer Style	1		1
Skin Staple Remover, Scissor Style		1	
Antiseptics			
PVP Prep Pad	1	1	
Quantity/Case	50	50	48



Suture Removal Trays

- Sharp instrument every time
- Fine, metal-tip scissors provide easy access to suture
- Convenient easy-to-open packaging
- Packaged in typical order of use



Suture Removal Trays	4366SUT1001	4366SUT1002	4366SUT1003	4366SUT1004	4366SUT1005	6422SUT1006	6422SUT1007	6422SUT1008	6422SUT1009	4366SUT1010	4366SUT1011
Packaging											
Hard Tray w/ Tyvek® Lid	1	1	1	1	1	1	1	1	1	1	1
Prep and Dressing											
Gauze, 3" x 3", 8-ply									1		
Gauze, 3" x 3", 12-ply	1	1	1	1	1	1	1	1		1	1
Instruments											
Scissors, Iris, Straight, 4½"	1	1		1						1	1
Scissors, Iris, Straight, 4½", Satin											1
Scissors, Littauer, 4½"					1						
Scissors, Littauer, 4½", Wireform						1	1				
Scissors, Littauer, 4¾", Plastic								1			
Forceps, Thumb, 4", Wireform		1					1		1		
Forcep, Adson, Serrated, 4¾"			1	1	1					1	
Forcep, Adson, Serrated, 4¾", Satin											1
Forceps, Finepoint, 4½", Plastic	1										
Forceps, Tissue, Plastic						1		1			
Antiseptics											
Alcohol Prep Pad, Large		1	1	1	1						
PVP Prep Pad		1	1	1	1						
Quantity/Case	50	50	50	50	50	50	50	50	50	50	50



Tracheostomy Trays

- Highly absorbent pre-split dressing
- Flexible trach brush for easy passage through inner cannula to minimize irritation
- Angle-cut twill ties to aid in easy threading through cannula slots
- Convenient easy-to-open packaging
- Packaged in typical order of use



Tracheostomy Trays	4366TRAC1001	4366TRAC1002	6422TRAC1003	6422TRAC1004	6422TRAC1005
Packaging					
Hard Tray, 3 Compartment w/Tyvek® Lid	1	1			
Hard Tray, w/ Tyvek® Lid			1	1	1
Prep and Dressing					
Gauze, Split, 4" x 4", 12-ply	4	4			
Gauze 4" x 4", 12-ply				3	4
Trach Sponge, 4" x 4"				1	1
Trach Dressing, Pre-Sewn	1	1			
Miscellaneous					
Cotton-Tip Applicator, 4", Plastic				2	
Cotton-Tip Applicator, 6", Plastic	3	3	3		2
Forcep, Plastic, Finepoint			1		
Pipe Cleaner	3	3	3	4	4
Twill Tape, 1/2" x 36"	1	1	1	1	1
Tracheostomy Brushes	1	1	1	1	1
Hydrogen Peroxide Packet, 1 oz.		2			
Gloves, LF/PF, Cuffed Medium, Pair	1	1	1	1	1
Drape, 13" x 19", 2-ply, T/P, Flat	1	1			
Drape, 18" x 15", Blue, 2/Pk			2	1	1
Quantity/Case	20	20	20	20	20



Fetal Monitoring Belts

- Soft woven polyester belts stretch easily
- Available in buttonhole with or without buttons, Velcro® style or mesh
- Buttonhole style is 2¾" by 48", with buttonholes every ½"
- Velcro® style is 1½" by 42", with 1" x 3" rounded ends to prevent nicks
- Packaged in pairs, one pink and one blue. Mesh belts are both white
- Disposable



Item No.	Description	Quantity
3596222362	Fetal Monitor Belt, White Mesh, 2 per pack	2 / PK, 100 / CA
3596223240	Fetal Monitor Belt, Velcro, 2 per pack	2 / PK, 100 / CA
3596224220	Fetal Monitor Belt, Button Holes, 2 per pack	2 / PK, 100 / CA
3596224221	Fetal Monitor Belt, Button Holes with Buttons, 2 per pack	2 / PK, 100 / CA



Baby Bonnets

- Knitted with 60% polyester, 40% cotton yarn for comfort and stretch
- 4" wide by 6" long
- Quality construction
- Hypo-allergenic
- Available in single-ply and double-ply



Item No.	Description	Quantity
3596099955	Baby Bonnet, Tri-Color, 2 Ply	1 / EA, 200 / CA
3596099956	Baby Bonnet, Tri-Color, 1 Ply	1 / EA, 200 / CA



Baby Lotion

- Gentle pH balanced formula with conditioners and emollients that hydrate sensitive skin
- Baby fresh scent
- Does not contain animal fat
- Convenient flip-top cap for easy dispensing
- Absorbs quickly and does not leave skin greasy
- Made in USA



Item No.	Description	Quantity
3306BL4020	Baby Lotion, 2oz.	1 / EA, 144 / CA
3306BL4040	Baby Lotion, 4oz.	1 / EA, 96 / CA

Complementary Products

Item Number	Description
3306PC8120	Tearless Baby Shampoo, 2 oz.



Baby Shampoo

- Tearless formula
- Rinses quickly and easily
- Leaves hair smelling fresh
- Gentle cleansing
- pH balanced
- Convenient flip-top cap for easy dispensing
- Made in USA



Item No.	Description	Quantity
3306PC8120	Tearless Baby Shampoo, 2oz., Flip Cap	1 / EA, 144 / CA



Perineal Cold Packs

- Extra-absorbent OB pad with adhesive strip to keep pack in place
- Supplies uniform therapeutic temperatures for post partum relief; helps reduce discomfort and swelling
- Designed for quick, easy activation for instant relief
- Degree Range: Cools to approximately 33°F / .5°C
- Maintains a safe therapeutic temperature for up to 1 hour
- Single-use



Item No.	Description	Quantity
1480069902	Perineal Cold Pack, 4" x 11"	1 / EA, 24 / CA



PERINEAL WASHES Perineal Wash Spray Bottle

- Individually bagged
- Graduated markings in both "oz" and "ml"
- Can be used for obstetric perineal care and ostomy care



Item No.	Description	Quantity
1233PPB1080	Perineal Irrigation Bottle, 8 oz., Bagged	1 / EA, 50 / CA



Biopsy Punches

- Seamless stainless steel blade with protective cap on punch blade
- Machined blade minimizes tissue tears
- High-impact polystyrene hollow handle for easy sample collection
- Punch size is printed on handle
- Package is color-coded by size for easy identification
- Single-use



Item No.	Description	Quantity
6422DP0200	Dermal Biopsy Punch, 2mm	1 / EA, 25 / BX
6422DP0300	Dermal Biopsy Punch, 3mm	1 / EA, 25 / BX
6422DP0350	Dermal Biopsy Punch, 3.5mm	1 / EA, 25 / BX
6422DP0400	Dermal Biopsy Punch, 4mm	1 / EA, 25 / BX
6422DP0500	Dermal Biopsy Punch, 5mm	1 / EA, 25 / BX
6422DP0600	Dermal Biopsy Punch, 6mm	1 / EA, 25 / BX
6422DP0800	Dermal Biopsy Punch, 8mm	1 / EA, 25 / BX



Disposable Safety Scalpels

- Highest quality stainless steel Sheffield blade, made in England
- Stable, fixed stainless steel blade with moveable sheath to prevent accidental blade exposure
- OSHA-compliant
- Protective, transparent sheath can be opened, closed and locked with a single hand; audible double-click when sheath is in locked position
- Sturdy, "easy-grip" acrylo-nitrile handle fits comfortably in either hand
- Sterile and individually packaged



Item No.	Description	Quantity
6422055610	Disposable Safety Scalpel #10, Sterile	1 / EA, 10 / BX
6422055611	Disposable Safety Scalpel #11, Sterile	1 / EA, 10 / BX
6422055615	Disposable Safety Scalpel #15, Sterile	1 / EA, 10 / BX



Esmark Bandages

- Smooth, comfortable material that ensures patient protection
- Consistent elongation for uniform compression and easy application
- Safely provides maximum pressure to extremities for exsanguination before surgery
- 6 sizes of both latex and latex-free varieties
- Conveniently packaged in sterile pouch



Item No.	Description	Quantity
6422055501	Esmark Bandage, Latex, 3" x 12', Sterile	1 / EA, 20 / CA
6422055502	Esmark Bandage, Latex, 3" x 9', Sterile	1 / EA, 20 / CA
6422055503	Esmark Bandage, Latex, 4" x 12', Sterile	1 / EA, 20 / CA
6422055504	Esmark Bandage, Latex, 4" x 9', Sterile	1 / EA, 20 / CA
6422055505	Esmark Bandage, Latex, 6" x 12', Sterile	1 / EA, 20 / CA
6422055506	Esmark Bandage, Latex, 6" x 9', Sterile	1 / EA, 20 / CA
6422055511	Esmark Bandage, Latex-Free, 3" x 12', Sterile	1 / EA, 20 / CA
6422055512	Esmark Bandage, Latex-Free, 3" x 9', Sterile	1 / EA, 20 / CA
6422055513	Esmark Bandage, Latex-Free, 4" x 12', Sterile	1 / EA, 20 / CA
6422055514	Esmark Bandage, Latex-Free, 4" x 9', Sterile	1 / EA, 20 / CA
6422055515	Esmark Bandage, Latex-Free, 6" x 12', Sterile	1 / EA, 20 / CA
6422055516	Esmark Bandage, Latex-Free, 6" x 9', Sterile	1 / EA, 20 / CA



SKIN MARKERS Correct Site Non-Sterile Skin Markers

- Assists with Joint Commission's "Universal Protocol for Preventing Wrong Site, Wrong Procedure, Wrong Person Surgery"
- Mini-size skin markers are designed to reduce cost and biowaste generation and encourage disposal after single patient use
- Full-size barrel diameter for hand comfort and better grip
- Mini skin markers contain Gentian violet ink and proof-of-use packaging
- Skin markers available in individual poly sleeve or dispenser box
- Temporary skin laminate tattoos cannot be moved or altered once placed (Joint Commission Compliant), yet can be removed without leaving marks
- Tattoos packaged in clear plastic box, 10 tattoos per poly zip bag



Item No.	Description	Quantity
4808MKR2100	Correct Site Skin Marker, Mini, Non Sterile, Dispenser Box of 100 Pens	1 / EA, 1200 / CS
4808MKR2250	Correct Site Skin Marker, Mini, Non Sterile w/Sleeve, 250/BX	1 / EA, 1250 / CS
4808MKR3001A	Correct Site Skin Tattoo, Small, YES, 1/2" Diameter, 100/BX, 4 BX/CS	1 / EA, 400 / CS
4808MKR3002A	Correct Site Skin Tattoo, Large, YES, 1" Diameter, 100/BX, 4 BX/CS	1 / EA, 400 / CS




Suction Tubing

- Durable and flexible tubing withstands heavy-duty suctioning
- Available in 6', 10', 12' and 20' lengths
- Soft female connector on each end
- Male adapter included
- Non-conductive, sterile PVC
- Suction power is 22" of mercury



Item No.	Description	Quantity
0205ST1001	Suction Tubing, 2 Connectors with Male Adapter, 3/16" x 6'	1 / EA, 50 / CA
0205ST1002	Suction Tubing, 2 Connectors with Male Adapter, 3/16" x 10'	1 / EA, 50 / CA
0205ST1003	Suction Tubing, 2 Connectors with Male Adapter, 3/16" x 12'	1 / EA, 20 / CA
0205ST1004	Suction Tubing, 2 Connectors with Male Adapter, 3/16" x 20'	1 / EA, 20 / CA
0205ST2001	Suction Tubing, 2 Connectors with Male Adapter, 1/4" x 6'	1 / EA, 50 / CA
0205ST2002	Suction Tubing, 2 Connectors with Male Adapter, 1/4" x 10'	1 / EA, 50 / CA
0205ST2003	Suction Tubing, 2 Connectors with Male Adapter, 1/4" x 12'	1 / EA, 20 / CA
0205ST2004	Suction Tubing, 2 Connectors with Male Adapter, 1/4" x 20'	1 / EA, 20 / CA
0205ST3001	Suction Tubing, 2 Connectors with Male Adapter, 9/32" x 6'	1 / EA, 50 / CA
0205ST3002	Suction Tubing, 2 Connectors with Male Adapter, 9/32" x 10'	1 / EA, 50 / CA
0205ST3003	Suction Tubing, 2 Connectors with Male Adapter, 9/32" x 12'	1 / EA, 20 / CA
0205ST3004	Suction Tubing, 2 Connectors with Male Adapter, 9/32" x 20'	1 / EA, 20 / CA



Tubing Connectors

- Sturdy PVC connectors provide tight seal
- Easy to connect to a variety of tubing sizes
- Available in sterile or non-sterile



Item No.	Description	Quantity
0205CON1001	Sims Tubing Connector, Sterile	1 / EA, 100 / CA
0205CON1002	Sims Tubing Connector, Non-Sterile	1 / EA, 100 / CA
0205CON2001	5-in-1 Tubing Connector, Sterile	1 / EA, 100 / CA
0205CON2002	5-in-1 Tubing Connector, Non-Sterile	1 / EA, 100 / CA



Yankauers

- PVC suction instrument
- Transparent plastic allows visual check for fluids
- Available with regular open tip or bulb tip with or without a control vent
- Curve in the Yankauer allows for ease of site use and minimizes hand and arm activity
- Non-conductive



Item No.	Description	Quantity
0205YANK1001	Yankauer, Bulb Tip, without Vent	1 / EA, 50 / CA
0205YANK1002	Yankauer, Bulb Tip, with Vent	1 / EA, 50 / CA
0205YANK2001	Yankauer, Open Tip, without Vent	1 / EA, 50 / CA
0205YANK2002	Yankauer, Open Tip, with Vent	1 / EA, 50 / CA
0205YANK3001	Yankauer, Bulb Tip, without Vent, Rigid	1 / EA, 50 / CA



Alcohol and Peroxide

- First-aid antiseptic
- Prevents the risk of infection in cuts, scrapes, and burns
- Alcohol may also be used as a cooling and refreshing rubbing/massaging compound
- Peroxide can be used as an oral debriding agent
- Peroxide also aids in the removal of phlegm and mucus associated with sore mouth



Item No.	Description	Quantity
3306HP0304	Hydrogen Peroxide 3% 4 oz.	1 / EA, 24 / CA
3306HP0308	Hydrogen Peroxide 3% 8 oz.	1 / EA, 12 / CA
3306HP0316	Hydrogen Peroxide 3% 16 oz.	1 / EA, 12 / CA
3306IA7016	Isopropyl Alcohol 70% 16 oz.	1 / EA, 12 / CA



Instant Hand Sanitizers

- Kills 99.9% of germs
- Enriched with aloe vera, Vitamin E and other emollients to condition and moisturize skin during repeated use
- Gel contains 62% ethyl alcohol
- Meets CDC standards on proper hand washing protocols
- Pleasant, honey fragrance
- Hypoallergenic
- Made in USA



Item No.	Description	Quantity
3306088300	Hand Sanitizer Alcohol Gel 4 oz.	1 / EA, 96 / CA
3306088303	Hand Sanitizer Alcohol Gel 8 oz. w/Pump	1 / EA, 48 / CA



Instant Cold Packs

- For temporary relief of minor pain and swelling caused by sprains, minor burns, contusions, or tooth aches and for use whenever cold therapy is indicated
- A simple squeeze and shake activates the pack
- Safety protocol recommends that the pack be wrapped in towel or wash cloth prior to use
- Available in 2 sizes—Large: 6" x 9", Junior: 5" x 7"
- Single patient use
- Disposable



Item No.	Description	Quantity
1480069903	Instant Cold Pack Junior 5" x 7"	1 / EA, 24 / CA
1480069904	Instant Cold Pack Large 6" x 9"	1 / EA, 24 / CA



Instant Warm Packs

- Soothing heat therapy for minor muscle aches, joint pain, stiffness and muscle spasms or whenever warm therapy is needed
- A simple squeeze and quick shake activates the pack
- Warms to safe 120°F
- Convenient single-use disposable compress
- DEHP-free
- Available in flock covered pack or laminated pack with protective sleeve (safety protocol recommends that poly pack be wrapped in towel or cloth sleeve provided prior to use)
- Fast, safe and reliable



Item No.	Description	Quantity
1480069901	Instant Warm Pack with Sleeve 6" x 9"	1 / EA, 24 / CA
148069951FC	Flock Covered Instant Warm Pack 6" x 9"	1 / EA, 24 / CA



Reusable Hot or Cold Gel Packs

- Activate in microwave for warm therapy or freezer for cold therapy
- Made of a soft-freeze gel that will remain flexible to mold around patient's body
- Non-toxic, non-irritating and non-caustic
- Primary chemical component: Super Absorbent Polymer (SAP)
- Packs require a cover and are not recommended to be applied directly to the skin
- Single patient use



Item No.	Description	Quantity
1480069905	Reusable Hot or Cold Gel Pack, 6" x 10"	1 / EA, 24 / CA



Adult Wipes

- Large, 9" x 13" disposable wipe for incontinence, bed bathing or general cleansing
- 45-gram Spunlace material
- Contains Aloe and Lanolin to promote skin health
- Alcohol-free, hypoallergenic
- Available in single-dispensing soft packs and tubs
- Fresh scent or fragrance-free



Item No.	Description	Quantity
33660713MC	Adult Wipe, Premium, Pre-moistened, Unscented Soft Pack w/ Resealable Lid 9" x 13", 50/PK	1 / EA, 12 / CA
33660813MC	Adult Wipe, Premium, Pre-Moistened, Scented, Soft Pack Refill 9" x 13", 50/BG	50 / PK, 600 / CA
33660913MC	Adult Wipe, Premium, Pre-Moistened, Scented, Hard Tub 9" x 13", 60/tub	60 / BX, 720 / CA



Washcloths

- Soft material is an economical alternative to cloth
- Quick absorption, strong even when wet
- Air-laid washcloths manufactured with cellulose
- Thermabond washcloth manufactured with rayon, polypropylene and viscous material
- First-quality material (not factory seconds)
- Available in 10" x 13" and 12" x 13" dimensions
- Single-use, non-sterile



Item No.	Description	Quantity
3596071013	Washcloth Airlaid 10"x13" Medium Weight	1 / EA, 500 / CA
3596071213	Washcloth Airlaid 12"x13" Medium Weight	1 / EA, 500 / CA
3596101305	Premium Washcloth 10"x13" Thermo-bond	1 / EA, 500 / CA



Odor Eliminators

- New clear formula—can be used on fabrics, skin and hair
- Odor elimination for general patient, necrotic tissue and incontinence odors
- Eliminates odors—doesn't just mask them
- Two formulas: Fresh Air Scented and Unscented
- Available in 2 oz. and 8 oz. pump spray bottles



Item No.	Description	Quantity
40670E2240	Biological Odor Eliminator, 2 oz. Pump Spray Scented	1 / EA, 24 / CA
40670E2241	Biological Odor Eliminator, 2 oz. Pump Spray Unscented	1 / EA, 24 / CA
40670E8120	Biological Odor Eliminator, 8 oz. Pump Spray Scented	1 / EA, 12 / CA
40670E8121	Biological Odor Eliminator, 8 oz. Pump Spray Unscented	1 / EA, 12 / CA



Oral Swabs

- Extra long 6" handle to avoid cross contamination
- Available in unflavored or fresh mint flavor
- Plastic stem won't absorb oral secretions
- Variety of packaging options: bulk, 10/pack and 20/pack
- Ridged foam head helps loosen debris and stimulate oral tissue
- Made in Canada



Item No.	Description	Quantity
4067TR0100	Oral Swab, Fresh Mint Flavored, Individually Wrapped, Bulk Packaged, 1000/CS	1 / EA, 1000 / CA
4067TR0200	Oral Swab, Fresh Mint Flavored, 20/PK, 25 PK/ BX; 2 BX/CS	20 / PK, 1000 / CA
4067TR0300	Oral Swab, Fresh Mint Flavored, Individually Wrapped, 10/ PK, 25PK/BX; 4 BX/CS	1 / EA, 1000 / CA
4067TT0100	Oral Swab, Unflavored, Individually Wrapped, Bulk Packaged, 1000/CS	1 / EA, 1000 / CA
4067TT0200	Oral Swab, Unflavored, 20/PK, 25 PK/ BX; 2 BX/ CS	20 / PK, 1000 / CA
4067TT0300	Oral Swab, Unflavored, Individually Wrapped, 10/PK, 25 PK/BX; 4 BX/CS	1 / EA, 1000 / CA



Tourniquets

- Easy to apply and remove with non-slip textured surface and minimal curl to stay on patient's arm
- Latex tourniquets are 100% natural with low powder content
- Latex-free tourniquets are made from a strong polyisoprene material with no latex proteins, resins or fatty acids
- Latex-free tourniquets come in powder-free blue or low-powder ivory
- New Orange ALERT tourniquet for easy visibility and patient safety
- High tensile and tear strength

Item No.	Description	Quantity
3308TRN1185	Latex-free Tourniquet, 1" x 18" Blue Roll, 25/RL, 20 RL/CS	1 / EA, 500 / CA
0565TRN1184	Latex-free Tourniquet, 1" x 18" Blue Flat Pack, 250/BX, 4 BX/CS	1 / EA, 1000 / CA
0565TRN1183	Latex-free Tourniquet, 1" x 18" Blue Rolled and Banded, 100/BG, 10 BG/CS	1 / EA, 1000 / CA
0565TRN1182	Latex-free Tourniquet, 1" x 18", Orange ALERT, Bulk Bag, 250/BG, 4 BG/CS	1 / EA, 1000 / CA



Suction Catheters

- Soft, strong PVC suction catheters
- Available in coiled or straight style
- Beveled whistle tip and offset eyes to promote maximum evacuation and help reduce mucosal trauma
- All catheters have control vent
- Adult sizes (12 Fr, 14 Fr, 16 Fr and 18 Fr) are 22" long and graduated for measuring proper depth insertion
- Pediatric sizes (6 Fr, 8 Fr and 10 Fr) are 16" long and graduated for measuring proper depth insertion, and also have Delee or rounded tip



Item No.	Description	Quantity
02050RSC06D	Suction Catheter, Delee Tip, 6Fr, Graduated, Straight	1 / EA, 50 / CA
02050RSC08D	Suction Catheter, Delee Tip, 8Fr, Graduated, Straight	1 / EA, 50 / CA
02050RSC10D	Suction Catheter, Delee Tip, 10Fr, Graduated Straight	1 / EA, 50 / CA
02050RSC12	Suction Catheter, Whistle Tip, 12Fr, Graduated, Straight	1 / EA, 50 / CA
02050RSC14	Suction Catheter, Whistle Tip, 14Fr, Graduated, Straight	1 / EA, 50 / CA
02050RSC16	Suction Catheter, Whistle Tip, 16Fr, Graduated, Straight	1 / EA, 50 / CA
02050RSC18	Suction Catheter, Whistle Tip, 18Fr, Graduated, Straight	1 / EA, 50 / CA
0205RSC06CD	Suction Catheter, Delee Tip, 6Fr, Graduated, Coiled	1 / EA, 50 / CA
0205RSC08CD	Suction Catheter, Delee Tip, 8Fr, Graduated, Coiled	1 / EA, 50 / CA
0205RSC10CD	Suction Catheter, Delee Tip, 10Fr, Graduated, Coiled	1 / EA, 50 / CA
0205RSC12C	Suction Catheter, Whistle Tip, 12Fr, Graduated, Coiled	1 / EA, 50 / CA
0205RSC14C	Suction Catheter, Whistle Tip, 14Fr, Graduated, Coiled	1 / EA, 50 / CA



Suction Catheter Kits

- Soft, strong PVC coiled respiratory catheters
- Beveled whistle tip and offset eyes promote maximum evacuation and help reduce mucosal trauma
- All catheters have control vent
- Adult sizes (12 Fr and 14 Fr) are 22" long and graduated for measuring proper depth insertion
- Pediatric sizes (6 Fr, 8 Fr and 10 Fr) are 16" long and graduated for measuring proper depth insertion, and also have Delee or rounded tip
- Kits include pre-cuffed powder-free vinyl gloves and pop-up solution container in a "walleled" jacket



Item No.	Description	Quantity
0205RSC06KD	Suction Catheter Kit, Delee Tip, 6Fr, Graduated, Sterile	1 / EA, 50 / CA
0205RSC08KD	Suction Catheter Kit, Delee Tip, 8Fr, Graduated, Sterile	1 / EA, 50 / CA
0205RSC10KD	Suction Catheter Kit, Delee Tip, 10Fr, Graduated, Sterile	1 / EA, 50 / CA
0205RSC12K	Suction Catheter Kit, Whistle Tip, 12Fr, Graduated, Sterile	1 / EA, 50 / CA
0205RSC14K	Suction Catheter Kit, Whistle Tip, 14Fr, Graduated, Sterile	1 / EA, 50 / CA



Suction Catheter Trays

- Soft, strong PVC coiled respiratory catheters
- Beveled whistle tip and offset eyes promote maximum evacuation and help reduce mucosal trauma
- Pediatric sizes (6 Fr, 8 Fr, and 10 Fr) have rounded Delee tip
- Control vent
- Adult size 14 Fr, 22" long and graduated for measuring proper depth insertion
- Kits include pre-cuffed powder-free vinyl gloves, rigid solution basin and 100ML container of saline solution



Item No.	Description	Quantity
0205RSC06TSD	Suction Catheter Tray with Saline, 6Fr	1 / EA, 36 / CA
0205RSC08TSD	Suction Catheter Tray with Saline, 8Fr	1 / EA, 36 / CA
0205RSC10TSD	Suction Catheter Tray with Saline, 10Fr	1 / EA, 36 / CA
0205RSC12TS	Suction Catheter Tray with Saline, 12Fr	1 / EA, 36 / CA
0205RSC14TS	Suction Catheter Tray with Saline, 14Fr	1 / EA, 36 / CA



Irrigation Syringes

- Easy to use sterile syringe for irrigation or enteral feeding
- Designed for one-handed operation and optimal control
- Tip protector and small tube adapter included
- Top-quality polypropylene
- Capacity: 60cc with easy-to-read graduations
- Available in thumb ring piston or bulb syringe



Item No.	Description	Quantity
0205BULB011P	Bulb Irrigation Syringe 60cc Latex-Free	1 / EA, 50 / CA
0205PIST015P	Piston Irrigation Syringe 60cc Latex-Free	1 / EA, 50 / CA



Pediatric Urine Collectors

- Versatile, unisex design provides a comfortable fit for infant or toddler
- PVC construction; 100% latex-free
- 5-oz. capacity with graduations to 160ml



Item No.	Description	Quantity
0205PUC300	Sterile Pediatric Urine Collector 160ml	1 / EA, 50 / BX



Foley Insertion Trays

- Sterile compact kit provides everything required to insert a Foley catheter with ease
- 1 pair vinyl gloves, 1 waterproof drape, 1 fenestrated drape, 1 package of catheter lubricant
- 1 pack of 3 PVP swab sticks to help prevent urinary tract infections
- Pre-filled syringe with 10cc or 30cc of sterile water
- One 1,000ml outer basin
- Catheter not included



Item No.	Description	Quantity
0205INTRAY80	Sterile Foley Insertion Tray 10cc Pre-filled Syringe	1 / EA, 20 / CA
0205INTRAY90	Sterile Foley Insertion Tray 30cc Pre-filled Syringe	1 / EA, 20 / CA



Urethral Catheterization Trays

- 14 Fr catheter—plastic or red rubber available
- 1 pair vinyl gloves, 1 waterproof drape, 1 fenestrated drape, 1 pack catheter lubricant, (1) 4-oz. specimen container
- 1 pack of 3 PVP swab sticks to help prevent urinary tract infections
- 1,000cc outer basin



Item No.	Description	Quantity
0205UCT870	Urethral Cath Tray Plastic Cath 14fr	1 / EA, 20 / CA
0205UCT871	Urethral Cath Tray Red Rubber 14fr	1 / EA, 20 / CA



Cotton Balls

- 100% Cotton
- Soft and absorbent
- Designed for a variety of uses throughout the hospital
- Packed in recloseable poly bags
- Nonsterile



Item No.	Description	Quantity
0565CTN500M	Medium Cotton Balls 500/BG	500 / BG, 4000 / CA
0565CTN200L	Large Cotton Balls 200/BG	200 / BG, 2000 / CA



Cohesive Bandages

- Self-adherent wrap for even, consistent compression in retaining dressings, protecting wounds, immobilizing injuries or securing devices to patients
- Bandage only sticks to itself and provides consistent compression to help minimize injured areas
- Fabric provides breathability and patient comfort
- Slightly tacky feel eliminates slipping or loosening
- Easy edge to quickly locate starting point of bandage
- Tears easily from the roll



Item No.	Description	Quantity
0565CBN1001	Cohesive Bandage Tan 1" x 5 yd.	1 / EA, 72 / CA
0565CBN1002	Cohesive Bandage Tan 2" x 5yd.	1 / EA, 72 / CA
0565CBN1003	Cohesive Bandage Tan 3" x 5yd.	1 / EA, 72 / CA
0565CBN1004	Cohesive Bandage Tan 4" x 5yd.	1 / EA, 72 / CA
0565CBN1006	Cohesive Bandage Tan 6" x 5yd.	1 / RL, 48 / CA
0565CBN115CP	Cohesive Bandage Color Pack 1½" x 5yd.	1 / EA, 72 / CA
0565CBN102CP	Cohesive Bandage Color Pack 2" x 5yd.	1 / EA, 72 / CA
0565CBN1103	Sterile Cohesive Bandage Tan 3" x 5yd.	1 / RL, 36 / CA
0565CBN1104	Sterile Cohesive Bandage Tan 4" x 5yd.	1 / EA, 36 / CA
0565CBN1106	Sterile Cohesive Bandage Tan 6" x 5yd.	1 / EA, 36 / CA



ELASTIC BANDAGES
Knit Elastic Bandages

- Cool and comfortable washable cotton bandage
- Knitted construction provides superior absorbency and breathability
- Heavy-gauge yarn and cotton-covered spandex for superior tensile strength and compression
- Resistant to deterioration from ointments and medications
- Four color-coded widths: 2", 3", 4", and 6"
- Available with dual self-closure flaps or metal clip closures



Item No.	Description	Quantity
3596REB2012	Standard Knit Elastic Bandage 2" x 5yd, Dual Self Closure 10/BX, 100/CA	1 / EA, 100 / CA
3596REB2013	Standard Knit Elastic Bandage 3" x 5yd, Dual Self Closure 10/BX, 100ea/ca	1 / EA, 100 / CA
3596REB2014	Standard Knit Elastic Bandage 4" x 5yd, Dual Self Closure 10/BX, 100ea/ca	1 / EA, 100 / CA
3596REB2016	Standard Knit Elastic Bandage 6" x 5yd, Dual Self Closure 10/BX, 100/CS	1 / EA, 100 / CA
3596REB2022	Standard Knit Elastic Bandage 2" x 5yd, Clip Closure 10/BX, 100/CS	1 / EA, 100 / CA
3596REB2023	Standard Knit Elastic Bandage 3" x 5yd, Clip Closure 10/BX, 100/CS	1 / EA, 100 / CA
3596REB2024	Standard Knit Elastic Bandage 4" x 5yd, Clip Closure 10/BX, 100/CS	1 / EA, 100 / CA
3596REB2026	Standard Knit Elastic Bandage 6" x 5yd, Clip Closure 10/BX, 100/CS	1 / EA, 100 / CA
3596REB3012	Premium Knit Elastic Bandage 2" x 5yd, Dual Self Closure 10/BX, 100/CS	1 / EA, 100 / CA

Item No.	Description	Quantity
3596REB3013	Premium Knit Elastic Bandage 3" x 5yd, Dual Self Closure 10/BX, 100/CS	1 / EA, 100 / CA
3596REB3014	Premium Knit Elastic Bandage 4" x 5½yd, Dual Self Closure 10/BX, 100/CS	1 / EA, 100 / CA
3596REB3016	Premium Knit Elastic Bandage 6" x 5½yd, Dual Self Closure 10/BX, 100/CS	1 / EA, 100 / CA
3596REB3112	Premium Knit Sterile Elastic Bandage 2" x 5yd, Dual Self Closure 20/CS	1 / EA, 20 / CA
3596REB3113	Premium Knit Sterile Elastic Bandage 3" x 5yd, Dual Self Closure 20/CS	1 / EA, 20 / CA
3596REB3114	Premium Knit Sterile Elastic Bandage 4" x 5½yd, Dual Self Closure 20/CS	1 / EA, 20 / CA
3596REB3116	Premium Knit Sterile Elastic Bandage 6" x 5½yd, Dual Self Closure 20/CS	1 / EA, 20 / CA



ELASTIC BANDAGES
Woven Elastic Bandages

- Cool and comfortable high cotton bandage absorbs more moisture than synthetic bands
- Economy, standard and premium bandages made of medium- to heavy-gauge yarns
- Orthopedic bandages made of heavy-gauge yarns and cotton-covered spandex and are washable
- Premium and orthopedic bandages resistant to deterioration from ointments and medications
- Metal clip closure
- Four widths: 2", 3", 4", and 6"



Item No.	Description	Quantity
3596REB1002	Economy Woven Elastic Bandage 2" x 4½yd, Clip Closure 10/BX, 100/CS	1 / EA, 100 / CA
3596REB1003	Economy Woven Elastic Bandage 3" x 4½yd, Clip Closure 10/BX, 100/CS	1 / EA, 100 / CA

contact your local Owens & Minor representative for details

Item No.	Description	Quantity
3596REB1004	Economy Woven Elastic Bandage 4" x 4½yd, Clip Closure 10/BX, 100/CS	1 / EA, 100 / CA
3596REB1006	Economy Woven Elastic Bandage 6" x 4½yd Clip Closure 10/BX, 100/CS	1 / EA, 100 / CA
3596REB2002	Standard Woven Elastic Bandage 2" x 5yd Clip Closure 10/BX, 100/CS	1 / EA, 100 / CA
3596REB2003	Standard Woven Elastic Bandage 3" x 5yd Clip Closure 10/BX, 100/CS	1 / EA, 100 / CA
3596REB2004	Standard Woven Elastic Bandage 4" x 5yd Clip Closure 10/BX, 100/CS	1 / EA, 100 / CA
3596REB2006	Standard Woven Elastic Bandage 6" x 5yd Clip Closure 10/BX, 100/CS	1 / EA, 100 / CA
3596REB2102	Standard Woven Sterile Elastic Bandage 2" x 5yd Clip Closure 20/CS	1 / EA, 20 / CA
3596REB2103	Standard Woven Sterile Elastic Bandage 3" x 5yd Clip Closure 20/CS	1 / EA, 20 / CA
3596REB2104	Standard Woven Sterile Elastic Bandage 4" x 5yd Clip Closure 20/CS	1 / EA, 20 / CA
3596REB2106	Standard Woven Sterile Elastic Bandage 6" x 5yd Clip Closure 20/CS	1 / EA, 20 / CA
3596REB3002	Premium Woven Elastic Bandage 2" x 5½yd Clip Closure 10/BX, 100/CS	1 / EA, 100 / CA
3596REB3003	Premium Woven Elastic Bandage 3" x 5½yd Clip Closure 10/BX, 100/CS	1 / EA, 100 / CA
3596REB3004	Premium Woven Elastic Bandage 4" x 5½yd Clip Closure 10/BX, 100/CS	1 / EA, 100 / CA
3596REB3006	Premium Woven Elastic Bandage 6" x 5½yd Clip Closure 10/BX, 100/CS	1 / EA, 100 / CA
3596REB3102	Premium Woven Sterile Elastic Bandage 2" x 5½yd Clip Closure 10/CS	1 / EA, 10 / CA
3596REB3103	Premium Woven Sterile Elastic Bandage 3" x 5½yd Clip Closure 10/CS	1 / EA, 10 / CA
3596REB3104	Premium Woven Sterile Elastic Bandage 4" x 5½yd Clip Closure 10/CS	1 / EA, 10 / CA
3596REB3106	Premium Woven Sterile Elastic Bandage 6" x 5½yd Clip Closure 10/CS	1 / EA, 10 / CA
3596REB4002	Orthopedic Woven Elastic Bandage 2" x 5yd Clip Closure 10/BX, 50/CS	1 / EA, 50 / CA
3596REB4003	Orthopedic Woven Elastic Bandage 3" x 5yd Clip Closure 10/BX, 50/CS	1 / EA, 50 / CA
3596REB4004	Orthopedic Woven Elastic Bandage 4" x 5yd Clip Closure 10/BX, 50/CS	1 / EA, 50 / CA
3596REB4006	Orthopedic Woven Elastic Bandage 6" x 5yd Clip Closure 10/BX, 50/CS	1 / EA, 50 / CA



Tubular Elastic Support Bandages

- Tubular elastic cotton/spandex bandage for support and compression
- Ideal for the effective management of sprains and strains, soft-tissue injuries, general edema, and post cast removal
- Easy to apply and stays in place without pins, clips or tape
- Not bulky which allows freedom of movement and comfort
- Delivers consistent support with a single application
- Convenient alternative to flat wrap bandages



Item No.	Description	Quantity
4350TSB5150	Tubular Elastic Support Bandage Natural Size A 1½" x 11yd, Infant Feet and Arms	1 / EA, 1 / CA
4350TSB5250	Tubular Elastic Support Bandage Natural Size B 2½" x 11yd, Small Hands and Arms	1 / EA, 1 / CA
4350TSB5275	Tubular Elastic Support Bandage Natural Size C 2¾" x 11yd, Medium Arms, Small Ankles	1 / EA, 1 / CA
4350TSB5300	Tubular Elastic Support Bandage Natural Size D 3" x 11yd, Large Arms, Medium Ankles, Small Knees	1 / EA, 1 / CA
4350TSB5350	Tubular Elastic Support Bandage Natural Size E 3½" x 11yd, Large Ankles, Medium Knees, Small Thighs	1 / EA, 1 / CA
4350TSB5400	Tubular Elastic Support Bandage Natural Size F 4" x 11yd, Large Knees, Medium Thighs	1 / EA, 1 / CA
4350TSB5450	Tubular Elastic Support Bandage Natural Size G 4½" x 11yd, Large Thighs	1 / EA, 1 / CA
4350TSB5675	Tubular Elastic Support Bandage Natural Size J 6¾" x 11yd, Small Trunks	1 / EA, 1 / CA
4350TSB5825	Tubular Elastic Support Bandage Natural Size K 8¾" x 11yd, Medium Trunks	1 / EA, 1 / CA

Item No.	Description	Quantity
4350TSB5013	Tubular Elastic Support Bandage Natural Size L 13" x 11yd, Large Trunks	1 / EA, 1 / CA
4350TSB7250	Tubular Elastic Support Bandage Beige Size B 2½" x 11yd, Small Hands and Arms	1 / EA, 1 / CA
4350TSB7275	Tubular Elastic Support Bandage Beige Size C 2¾" x 11yd Medium Arms, Small Ankles	1 / EA, 1 / CA
4350TSB7300	Tubular Elastic Support Bandage Beige Size D 3" x 11yd Large Arms, Medium Ankles, Small Knees	1 / EA, 1 / CA
4350TSB7350	Tubular Elastic Support Bandage Beige Size E 3½" x 11yd Large Ankles, Medium Knees, Small Thighs	1 / EA, 1 / CA
4350TSB7400	Tubular Elastic Support Bandage Beige Size F 4" x 11yd Large Knees, Medium Thighs	1 / EA, 1 / CA
4350TSB7450	Tubular Elastic Support Bandage Beige Size G 4½" x 11yd, Large Thighs	1 / EA, 1 / CA
4350TSB9250	Single Unit Tubular Elastic Support Bandage Natural Size B 2½" x 36" Small Hands and Arms	1 / EA, 12 / CA
4350TSB9275	Single Unit Tubular Elastic Support Bandage Natural Size C 2¾" x 36" Medium Arms, Small Ankles	1 / EA, 12 / CA
4350TSB9300	Single Unit Tubular Elastic Support Bandage Natural Size D 3" x 36" Large Arms, Medium Ankles, Small Knees	1 / EA, 12 / CA
4350TSB9350	Single Unit Tubular Elastic Support Bandage Natural Size E 3½" x 36" Large Ankles, Medium Knees, Small Thighs	1 / EA, 12 / CA
4350TSB9400	Single Unit Tubular Elastic Support Bandage Natural Size F 4" x 36" Large Knees, Medium Thighs	1 / EA, 12 / CA



Tubular Net Bandages

- Retainer net holds wound dressings, hot and cold packs, and medical devices in place without tape
- Open-weave design allows air flow and permits freedom of movement
- Cuts easily and stretches to fit various body parts and shapes
- 100% latex-free elastic and polyester
- Whiter and softer than competitive products
- Ravel-free



Item No.	Description	Quantity
4350NETBND01	Tubular Net Bandage Size 1 Fits Fingers, Toes, or Wrists	1 / RL, 10 / CA
4350NETBND02	Tubular Net Bandage Size 2 Fits Small Hand, Arm, Leg, or Foot	1 / RL, 10 / CA
4350NETBND03	Tubular Net Bandage Size 3 Fits Medium Hand, Arm, Leg, or Foot	1 / RL, 10 / CA
4350NETBND04	Tubular Net Bandage Size 4 Fits Large Hand, Arm, Leg, or Foot	1 / RL, 10 / CA
4350NETBND05	Tubular Net Bandage Size 5 Fits Small Head, Shoulder, or Thigh	1 / EA, 10 / CA
4350NETBND06	Tubular Net Bandage Size 6 Fits Large Head, Shoulder, or Thigh	1 / RL, 10 / CA
4350NETBND07	Tubular Net Bandage Size 7 Fits Small Chest, Back, Perineum, or Axilla	1 / RL, 10 / CA
4350NETBND08	Tubular Net Bandage Size 8 Fits Medium Chest, Back, Perineum, or Axilla	1 / RL, 10 / CA
4350NETBND09	Tubular Net Bandage Size 9 Fits Large Chest, Back, Perineum, or Axilla	1 / RL, 10 / CA
4350NETBND10	Tubular Net Bandage Size 10 Fits X-Large Chest, Back, Perineum, or Axilla	1 / RL, 10 / CA
4350NETBND11	Tubular Net Bandage Size 11 Fits XX-Large Chest, Back, Perineum, or Axilla	1 / RL, 10 / CA
4350NETBND12	Tubular Net Bandage Size 12 Fits 3X-Large Chest, Back, Perineum, or Axilla	1 / RL, 10 / CA
4350NETBND22	Tubular Net Bandage Size 22 Custom Size Stretches to 80"	1 / RL, 10 / CA



Paste Bandages

- 100% cotton gauze bandage used in the compression treatment of ulcers, edema and strains of the extremities
- Non-raveling gauze evenly impregnated with a non-hardening zinc oxide paste
- No gelatin or preservatives—no expiration date
- Bandages conform to patient's leg without gaps or wrinkles
- Available with calamine lotion to ease itching
- Width: 3" or 4"



Item No.	Description	Quantity
3596099903	Paste Bandage 3" x 10yd	1 / EA, 12 / CA
3596099903C	Paste Bandage 3" x 10yd with Calamine	1 / EA, 12 / CA
3596099904	Paste Bandage 4" x 10yd	1 / EA, 12 / CA
3596099904C	Paste Bandage 4" x 10yd with Calamine	1 / EA, 12 / CA



Compression Bandage System

- 4-Layer system for optimum compression and absorbency in the treatment of venous leg ulcers and related conditions
- Layer 1: Non-woven padding for comfort and protection of bony joints
- Layer 2: Crepe bandage which anchors and smoothes padding layer while absorbing exudate
- Layer 3: Light compression bandage that conforms to any contour while providing half of the compression of the system
- Layer 4: Cohesive bandage that sticks to itself and secures the system together
- Adhesive tape strips and application guide included



Item No.	Description	Quantity
3596BND9001	4-Layer Compression Bandage System Layers Individually Packed	1 / PK, 8 / CA

A

Adult and Pediatric Professional Stethoscopes 29
 Adult Wipes 55
 Airways, Berman 14
 Alcohol..... 53
 Aneroid Sphygmomanometers, Deluxe 26
 Aneroid Sphygmomanometers, Latex-Free Deluxe 26
 Aneroid Sphygmomanometers, Latex-Free Economy 26
 Aneroid Sphygmomanometers, Economy 25
 Antiperspirant Deodorant, Roll On 11

B

Baby Bonnets 46
 Baby Lotion..... 47
 Baby Shampoo 47
 Bandage System, Compression 66
 Bandages, Cohesive 62
 Bandages, Esmark 50
 Bandages, Knit Elastic 63
 Bandages, Paste 66
 Bandages, Tubular Elastic Support 64
 Bandages, Tubular Net 65
 Bandages, Woven Elastic 63
 Bed Pillows, Disposable 12
 Bed Pillows, Limited Reusable 12
 Bed Pillows, Reusable 13
 Belts, Fetal Monitoring..... 46
 Berman Airways 14
 Biopsy Punches..... 49
 Blanket Warmers, Counter-Top..... 33
 Blanket Warmers, Double-Door 34
 Blanket Warmers, Single-Door 34
 Blanket Warmers, Under-Counter 33
 Body Bath and Shampoo 11
 Body Bath and Shampoo, No Rinse 11
 Bonnets, Baby 46
 Briefs, Mesh 15

C

Cardiology, Dual-Head Stethoscopes..... 30
 Catheter, Suction Kits..... 58
 Catheter, Suction Trays 59
 Catheterization Trays, Urethral..... 61
 Catheters, Suction 58
 Central Line Dressing Trays 36
 Chart Paper, EKG 32
 Chrome and Plastic Disposable Stethoscopes 27
 Cleaner, Enzyme Gel Spray 23
 Cleaner, Enzyme Pre-Soak 23
 Cleaners, Neutral pH Instrument..... 23

Clear Polyethylene Telemetry Pouches 17
 Coats, Premium Disposable Lab 22
 Coats, Standard Disposable Lab..... 22
 Cohesive Bandages..... 62
 Cold Packs, Instant 54
 Cold Packs, Perineal 48
 Cold Reusable Gel Packs 55
 Collectors, Pediatric Urine 60
 Compression Bandage System..... 66
 Connectors, Tubing..... 51
 Correct Site Non-Sterile Skin Markers..... 50
 Cotton Balls..... 62
 Counter-Top Blanket Warmers..... 33
 Cover Gowns 18
 Coveralls, Premium Disposable 21
 Coveralls, Standard Disposable 21
 Cups, Denture 10

D

Deluxe Aneroid Sphygmomanometers..... 26
 Deluxe Aneroid Sphygmomanometers, Latex-Free 26
 Denture Cups 10
 Deodorant, Roll-On..... 11
 Dermal Biopsy Punches 49
 Digital Thermometers, Oral and Rectal 31
 Digital Thermometers, Oral Kits with Covers 31
 Disposable Chrome and Plastic Stethoscopes 27
 Disposable Lab Coats, Premium 22
 Disposable Lab Coats, Standard..... 22
 Disposable Pillows 12
 Disposable Plastic Stethoscopes..... 27
 Disposable Premium Coveralls 21
 Disposable Procedure Masks 18
 Disposable Safety Scalpels..... 49
 Disposable Standard Coveralls 21
 Disposable Washcloths 56
 Double-Door Blanket Warmers..... 34
 Drainage and Incision Trays 36
 Dressing Trays, Central Line 36
 Dual-Head Cardiology Stethoscopes..... 30
 Dual-Head Neonatal Stethoscopes 27
 Dual-Head Pediatric Stethoscopes..... 28
 Dual-Head, Infant Stethoscopes 28
 Dual-Head, Lightweight Stethoscopes..... 29

E

ECG Paper 32
 Economy Aneroid Sphygmomanometers..... 25
 Economy Aneroid Sphygmomanometers, Latex-Free 26
 EKG Paper 32
 Elastic Bandages, Woven 63

Elastic Bandages, Knit	63
Elastic Support Bandages, Tubular	64
Enzymatic Cleaner, Pre-Soak	23
Enzymatic Pre-Cleaner, Gel Spray	23
Enzyme Cleaner Pre-Soak	23
Enzyme Gel Spray Pre-Cleaner	23
Esmark Bandages	50

F

Fall Management Slippers	15
Fetal Monitoring Belts	46
Foam Telemetry Pouches	17
Foley Insertion Trays	61
Forcep, Thumb	24

G

Gel Packs, Reusable Hot or Cold	55
Gowns, Cover	18
Gowns, Level 1 Isolation	19
Gowns, Level 2 Isolation	20
Gowns, Level 3 Isolation	20
Gowns, Personal Protective	19

H

Hand and Body Lotion	10
Hand Sanitizers, Instant	53
Hot Reusable Gel Packs	55

I

Incision and Drainage Trays	36
Infant Dual-Head Stethoscopes	28
Insertion Trays, Foley	61
Instant Cold Packs	54
Instant Hand Sanitizers	53
Instant Warm Packs	54
Instrument Cleaners, Neutral pH	23
Instrument Lubricants	23
Irrigation Syringes	60
Irrigation Trays with Bulb Syringe	35
Irrigation Trays with Piston Syringe	35
Isolation Cover Gowns, Yellow	18
Isolation Gowns, Level 1	19
Isolation Gowns, Level 2	20
Isolation Gowns, Level 3	20
IV Start Kits	38
IV Start Kits with ChloraPrep	40
IV Start Kits with Chlorascrub	39

J

Jelly, Lubricating	25
--------------------------	----

K

Kits, IV Start	38
Kits, IV Start with ChloraPrep	40
Kits, IV Start with Chlorascrub	39
Kits, Suction Catheter	58
Knit Elastic Bandages	63

L

Lab Coats, Premium Disposable	22
Lab Coats, Standard Disposable	22
Latex-Free Deluxe Aneroid Sphygmomanometers	26
Latex-Free Economy Aneroid Sphygmomanometers	26
Level 1 Isolation Gowns	19
Level 2 Isolation Gowns	20
Level 3 Isolation Gowns	20
Lightweight Dual-Head Stethoscopes	29
Lightweight Single-Head Bowles Stethoscopes	29
Lightweight Single-Head Stethoscopes	28
Limited Reusable Pillows	12
Littauer Scissors	24
Lotion, Baby	47
Lotion, Hand and Body	10
Lubricants, Instrument	23
Lubricating Jelly	25

M

Markers, Correct Site Non-Sterile Skin	50
Masks, Procedure	18
Mesh Briefs	15
Monitoring Belts, Fetal	46

N

Neonatal Dual-Head Stethoscopes	27
Net Bandages, Tubular	65
Neutral pH Instrument Cleaners	23
Non-Woven Telemetry Pouches with Neck and Waist Straps	17
Non-Woven Telemetry Pouches with Neck Straps	17
No-Rinse Shampoo and Body Bath	11

O

Odor Eliminators	56
Oral and Rectal Digital Thermometers	31
Oral Swabs	57
Oral Thermometer Kits with Covers	31

P

Packs, Instant Cold 54
 Packs, Instant Warm 54
 Packs, Perineal Cold 48
 Packs, Reusable Hot or Cold..... 55
 Panties, Mesh..... 15
 Paper, EKG..... 32
 Paste Bandages 66
 Pediatric and Adult Professional Stethoscopes 29
 Pediatric Dual-Head Stethoscopes..... 28
 Pediatric Urine Collectors 60
 Perineal Cold Packs 48
 Perineal Wash Spray Bottle 48
 Peroxide..... 53
 Personal Protective Gowns..... 19
 Pillows, Disposable 12
 Pillows, Limited Reusable..... 12
 Pillows, Reusable 13
 Plastic and Chrome Disposable Stethoscopes 27
 Plastic, Disposable Stethoscopes 27
 Pouches, Non-Woven Telemetry with Neck and Waist Straps 17
 Pouches, Non-Woven Telemetry with Neck Straps..... 17
 Pouches, Telemetry Clear Polyethylene 17
 Pouches, Telemetry Foam..... 17
 Pre-Cleaner, Enzyme Gel Spray 23
 Premium Disposable Coveralls 21
 Premium Disposable Lab Coats 22
 Pre-Soak, Enzyme Cleaner 23
 Procedure Masks 18
 Professional Adult and Pediatric Stethoscopes 29
 Protective Gowns, Personal..... 19
 Punches, Biopsy..... 49

Q

R

Rectal and Oral Digital Thermometers 31
 Removal Trays, Suture 44
 Remover Trays, Skin Staple 43
 Reusable Hot or Cold Gel Packs 55
 Reusable Pillows 13
 Roll-on Deodorant 11

S

Safety Scalpels, Disposable..... 49
 Sanitizers, Instant Hand 53
 Scalpels, Disposable Safety..... 49
 Scissors, Littauer 24
 Shampoo and Body Bath 11

Shampoo and Body Bath, No-Rinse 11
 Shampoo, Baby 47
 Sheaths, Thermometer 31
 Single-Door Blanket Warmers 34
 Single-Head, Bowles Lightweight Stethoscopes..... 29
 Single-Head, Lightweight Stethoscopes 28
 Skin Markers, Correct Site Non-Sterile..... 50
 Skin Staple Remover Trays 43
 Slippers, Fall Management 15
 Slippers, Smooth Side Out 16
 Slippers, Textured Side Out..... 16
 Smooth Side Out Slippers..... 16
 Sphygmomanometers, Deluxe Aneroid 26
 Sphygmomanometers, Economy Aneroid..... 25
 Sphygmomanometers, Latex-Free Deluxe Aneroid..... 26
 Sphygmomanometers, Latex-Free Economy Aneroid 26
 Sprague-Rappaport Stethoscopes with Accessory Kit..... 30
 Spray, Perineal Wash Bottle 48
 Standard Disposable Coveralls 21
 Standard Disposable Lab Coats..... 22
 Stethoscopes, Disposable Chrome and Plastic 27
 Stethoscopes, Dual-Head Cardiology..... 30
 Stethoscopes, Infant Dual-Head..... 28
 Stethoscopes, Lightweight Dual-Head 29
 Stethoscopes, Lightweight Single-Head 28
 Stethoscopes, Lightweight Single-Head Bowles 29
 Stethoscopes, Neonatal Dual-Head 27
 Stethoscopes, Pediatric Dual-Head 28
 Stethoscopes, Plastic Disposable 27
 Stethoscopes, Professional Adult and Pediatric 29
 Stethoscopes, Sprague-Rappaport with Accessory Kit..... 30
 Suction Catheter Kits..... 58
 Suction Catheter Trays 59
 Suction Catheters 58
 Suction Tubing..... 51
 Suction Yankauers 52
 Support Bandages, Tubular Elastic..... 64
 Suture Removal Trays 44
 Suture Trays 42
 Swabs, Oral 57
 Syringes, Irrigation..... 60

T

Telemetry Pouches, Clear Polyethylene 17
 Telemetry Pouches, Foam..... 17
 Telemetry Pouches, Non-Woven with Neck and Waist Straps 17
 Telemetry Pouches, Non-woven with Neck Straps 17
 Textured Side Out Slippers..... 16
 Thermometers, Oral and Rectal Digital..... 31
 Thermometers, Oral Kits with Covers 31
 Thermometers, Sheaths 31
 Thumb Forcep 24

Tourniquets	57
Tracheostomy Trays.....	45
Trays, Central Line Dressing	36
Trays, Foley Insertion	61
Trays, Incision and Drainage.....	36
Trays, Irrigation with Bulb Syringe.....	35
Trays, Irrigation with Piston Syringe	35
Trays, Skin Staple Remover	43
Trays, Suction Catheter	59
Trays, Suture.....	42
Trays, Suture Removal	44
Trays, Tracheostomy.....	45
Trays, Urethral Catheterization	61
Tubing Connectors.....	51
Tubular Elastic Support Bandages.....	64
Tubular Net Bandages.....	65

U

Under-Counter Blanket Warmers	33
Unna Boot.....	66
Urethral Catheterization Trays.....	61
Urine Collectors, Pediatric	60

V

W

Warm Packs, Instant	54
Warmers, Counter-Top Blanket.....	33
Warmers, Double-Door Blanket.....	34
Warmers, Single-Door Blanket.....	34
Warmers, Under-Counter Blanket.....	33
Wash, Perineal Spray Bottle	48
Washcloths.....	56
Wipes, Adult	55
Woven Elastic Bandages	63

X

Y

Yankauers.....	52
Yellow Cover Gowns.....	18
Yellow Isolation Gowns, Level 1.....	19
Yellow Isolation Gowns, Level 2.....	20
Yellow Isolation Gowns, Level 3.....	20

Z

The O&M Distribution Network

MEDICHOICE®
from Owens & Minor

HOME OFFICE

9120 Lockwood Boulevard
Mechanicsville, VA 23116
804.723.7000 • f 804.723.7100

ALBUQUERQUE

8000 Bluewater Road NW, Suite A
Albuquerque, NM 87121
505.352.1300 • f 505.352.3800

ALLENTOWN, PA

7437 Industrial Boulevard
Allentown, PA 18106-9378
610.706.3149 • f 610.391.9480

ATLANTA

3375 Cobb International Blvd.
Kennesaw, GA 30152-4399
770.426.3700 • f 770.426.3730

AUGUSTA

777 Horizon South Pkwy.
Grovetown, GA 30813-9043
706.860.3313 • f 706.860.2729

BALTIMORE / WASHINGTON

1425 Magellan Road
Hanover, MD 21076
410.865.7700 • f 410.691.2364

BIRMINGHAM

289 Cahaba Valley Parkway
Pelham, AL 35124-1146
205.403.4900 • f 205.403.4908

BOSTON

135 Constitution Boulevard
Franklin, MA 02038-2584
508.520.2000 • f 508.520.2080

CHARLOTTE

202 Commerce Boulevard
Kings Mountain, NC 28086
704.730.8117 • f 704.730.1514

CHICAGO | HANOVER PARK

945 Muirfield Drive
Hanover Park, IL 60133
847.273.1003 • f 630.307.0395

CHICAGO | TINLEY PARK

Hickory Creek Corporate Ctr. #1
7979 West 183rd Street
Tinley Park, IL 60477
708.429.8880 • f 708.429.5845

CINCINNATI

1605 Worldwide Boulevard
Hebron, KY 41048
859.282.7910 • f 859.282.7912

CLEVELAND

7905 Cochran Rd., Suite 300
Glenwillow, OH 44139-5471
440.715.1146 • f 440.914.9900

DALLAS

550 Lakeside Parkway, Suite 200
Flower Mound, TX 75028
972.538.5500 • f 972.393.7627

DENVER

11440 E. 56th Avenue
Denver, CO 80239-2228
303.576.4200 • f 303.375.1859

DES MOINES

4488 112th Street
Urbandale, IA 50322
515.276.0206 • f 515.251.8367

DETROIT

41133 Van Born Road, Suite 100
Belleville, MI 48111-1147
734.394.0575 • f 734.394.3037

FT. LAUDERDALE

14599 NW 8th Street
Sunrise, FL 33325-6248
954.845.2200 • f 954.845.2250

GAINESVILLE, FL

4807 NE 63rd Avenue
Gainesville, FL 32609
352.264.9939 • f 352.264.2289

GRAND RAPIDS, MI

3900 S. Greenbrooke Drive, SE
Kentwood, MI 49512-5326
616.698.7717 • f 616.698.7843

GREENSBURG, PA

134 Industrial Park Road
Suite 1000
Greensburg, PA 15601-0339
724.837.5700 • f 724.830.7025

HARLINGEN, TX

Valley Industrial Center #1
1805 N. Loop 499, Suite 130
Harlingen, TX 78550
956.428.6401 • f 956.428.8857

HONOLULU

91-240 Kalaeloa Boulevard
Kapolei, HI 96707
808.682.2444 • f 808.682.0555

HOUSTON

2700 Brittmoore
Houston, TX 77043-1402
713.939.1464 • f 713.939.3325

INDIANAPOLIS

6585 East 30th St., Suite A
Indianapolis, IN 46219-1101
317.568.4500 • f 317.549.1801

JACKSON

326 U.S. Hwy. 49 South
Richland, MS 39218-9439
601.933.6280 • f 601.933.6292

JACKSONVILLE

8489 Westside Industrial Drive
Jacksonville, FL 32219
904.596.4600 • f 904.783.8376

KANSAS CITY

1800 Southern Road
Kansas City, MO 64120-1144
816.483.6613 • f 816.231.5902

KNOXVILLE

3551 Workman Road
Knoxville, TN 37921
865.588.0676 • f 865.588.0705

LAS VEGAS

2970 North Lamb Boulevard
Suite 105
Las Vegas, NV 89115
702.531.7464 • f 702.531.7469

LATHROP, CA

17720 Shideler Parkway, Suite B
Lathrop, CA 95330-8336
209.858.4001 • f 209.858.4020

LOS ANGELES

455 South Brea Canyon Road
City of Industry, CA 91789-3058
909.444.6500 • f 909.468.9770

LOUISVILLE

5000 Commerce Crossing Street
Suite 100
Louisville, KY 40229
502.810.1400 • f 502.961.6947

MEMPHIS

5295 Logistics Drive, Suite 101
Memphis, TN 38118
901.794.9448 • f 901.360.1222

MINNEAPOLIS

2151 County Road H2W
Mounds View, MN 55112-4729
763.785.9100 • f 763.785.2615

NASHVILLE

3308 Briley Park Blvd. South
Nashville, TN 37207
615.650.9200 • f 615.650.9311

NEW ORLEANS

310 James Drive East, Suite 100
Saint Rose, LA 70087
504.472.5783 • f 504.472.5750

OKLAHOMA CITY

6100 NW 2nd St., Suite 100
Oklahoma City, OK 73127-6523
405.495.1900 • f 405.789.7620

OMAHA

10288 South 134th Street
Omaha, NE 68138
402.455.4444 • f 402.896.0884

ORLANDO

301 Gills Drive, Suite 100
Orlando, FL 32824
407.299.5490 • f 407.541.5359

PHOENIX

8313 W. Pierce Street, Suite 100
Tolleson, AZ 85353
602.269.7121 • f 602.278.7883

PORTLAND

27695 SW 95th Avenue
Wilsonville, OR 97070-8215
503.685.9800 • f 503.682.0228

RALEIGH

4200 Global Street
Raleigh, NC 27610
919.231.6367 • f 919.231.6352

RICHMOND

11900 North Lakeridge Parkway
Ashland, VA 23005
804.752.5200 • f 804.798.6728

ROCHESTER

45 Vantage Point Drive
Rochester, NY 14624
585.352.7060 • f 585.352.7080

SALT LAKE CITY

2475 South 3200 West, Suite A
West Valley City, UT 84119-1236
801.972.8188 • f 801.973.0734

SAN DIEGO

4000 Ruffin Road, Suite H
San Diego, CA 92123-1849
858.560.4651 • f 858.560.4298

SAN FRANCISCO

6150 Las Positas Road
Livermore, CA 94551-7809
925.606.1133 • f 925.606.1401

SEATTLE

4301 West Valley Highway East
Suite 300
Sumner, WA 98390
253.372.1100 • f 253.852.6774

ST. LOUIS

15 Gateway Commerce Center Dr. W.
Edwardsville, IL 62025
618.931.4396 • f 618.931.4661

TULSA

8831 East Pine, Suite 100
Tulsa, OK 74115-5801
918.838.8012 • f 918.838.2688

WAUNAKEE, WI

201 Raemisch Road
Waunakee, WI 53597-9664
608.849.2600 • f 608.849.2603

WEST COAST SALES OFFICE

452-4 Sespe Avenue
Fillmore, CA 93015
805.524.0243 • f 805.524.0220

WEST DEPTFORD, NJ

1220 Forest Parkway
West Deptford, NJ 08066
856.423.9900 • f 856.423.7387

About MediChoice®

MediChoice® is Owens & Minor's brand of medical-surgical products and equipment that is growing to meet the changing needs of healthcare. The MediChoice® product portfolio reflects Owens & Minor's longstanding commitment to building strong and diverse supply chain relationships. Products we source from women-, minority- and veteran-owned businesses make up over one-third of the MediChoice® product line.

About Owens & Minor

Owens & Minor, Inc., a FORTUNE 500 company headquartered in Richmond, Virginia, is the leading distributor of national name-brand medical and surgical supplies and a healthcare supply chain management company. With a diverse product and service offering and distribution centers throughout the United States, the company serves hospitals, integrated healthcare systems, alternate care locations, group purchasing organizations, and the federal government.

Contact Information

To learn more about MediChoice® products contact your local Owens & Minor representative or visit us online at owens-minor.com.